

PROUD TO BE CANADIAN

BUY LOCAL SUPPORT LOCAL



**PENGUIN
HEALTH**

Enhancing Patient Care

Q2 PROMOTION 2025

CALL IN TO ORDER TODAY
905-901-1579

PROMO EXPIRES
JUNE 30, 2025

LUNA WHITENING TAKE HOME KIT

5+5

VISIT PAGE 14

SciCan OPTIM WIPES 33TB PH-SCI33W12

BUY 24



BUY 12

VISIT PAGE 5


ATLAS RAPID INCUBATOR

SAVE
100%

VISIT PAGE 1

ATLAS SHIELD TYPE 5 CHEMICAL INTEGRATORS AT-4321

5+5



VISIT PAGE 3

HARVARD MULTICHROME COMPULES 7100002

4+1



VISIT PAGE 20

AMW AIRPEAK HIGH-SPEED HANDPIECES CODES VARY

2+1



VISIT PAGE 26

WWW.PENGUINHEALTH.CO

905-901-1579 | ORDERS@PENGUINHEALTHGROUP.COM

Table of Contents

Penguin Sterilization	1 - 5
Penguin Disposables	6-11
Preventative Dentistry	12 - 14
Pediatric Crowns	15 - 16
Restorative Dentistry	17 - 31
Oral Surgical Dentistry	32
PPE	33
Capital Equipment	34- 39
Construction Promo	

JOIN OUR PENGUIN+ REWARDS PROGRAM

PER ORDER

- ⊕ Spend \$1,000 per order Get a \$50 Gift Card
- ⊕ Spend \$2,000 per order Get a \$150 Gift Card
- ⊕ Spend \$3,500 per order Get a FREE B-Cure Curing Light! (Worth \$699)
- ⊕ Spend \$5,000 per order Get 2 FREE AMW Slow-Speed Handpieces (Worth \$800) or 1 FREE AMW High Speed Handpiece! (Worth \$800)

ANNUAL

- ⊕ Spend \$10,000 Get a FREE BEYES CanaView IntraOral Camera (WORTH \$1199)
- ⊕ Spend \$20,000 Get a FREE PENGUIN i-Sensor (WORTH \$4499!)
- ⊕ Spend \$30,000 Get a FREE Dental Wall-Mounted X-Ray (Worth \$5,500!)
- ⊕ Spend \$40,000 Get a FREE PT-A Air Polisher & Ultrasonic Scaler! (Worth \$8,000!)
- ⊕ Spend \$60,000 Get a FREE Dental Chair Package (Worth \$15,000+! - Installation & Freight Separate)

PROGRAM RULES



Can not be combined with other existing offers. Spending must be on products from our catalogues or promotions. Reward is distributed once target is reached. Program effective date July 1st, 2024 [Can not be back-dated]. Any form of special-pricing is to the discretion of the representative and may disqualify. Must choose either per-order or annual rewards when processing an order. Taxes, shipping and handling costs cannot be part of the total spend. All other company terms and conditions apply.



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.co




A PROUD CANADIAN
COMPANY

ATLAS RAPID INCUBATOR

SPORE TEST RESULTS IN
20 MINUTES!

SET UP YOUR OFFICE ON
ATLAS RAPID TODAY!

- Automate your data - no more human error
- Ultra rapid results in only 20-MINUTES! - no more waiting around and extensive quarantines



~~\$1499~~

WITH BUY BACK PROGRAM

SAVE 100%

#REORDER CODE **PH01980**



NET \$295.99/BX

4+1

REG \$369.99/BX

#REORDER CODE **AX224**

BUY BACK PROGRAM

Trade-In Your Old Incubator for a credit towards the new unit!
Save up to 100%

- Trade-in your OLD Incubator and Left-Over Spore tests
- Get our brand new \$1499-worth ATLAS Rapid Incubator for free

















CLICK OR SCAN to Watch the Video!



CLICK OR SCAN this QR Code for More Information

WHY MAKE THE SWITCH TO ATLAS RAPID? See the difference!

Benefits	24- Hour System	ATLAS Rapid 20-Minutes System
 Test Duration	24 hours	20 minutes
 Result Receipt	Manual	Automated Printout, minimizing human error
 Ease of Use	Varies	Touch screen display for simple operation and training
 Contamination Control	Varies	Lid-cover for spore tests to reduce contamination
 Cost	Standard	30% less expensive than leading brands
 Efficiency	Lower	Increased efficiency, more emergency patient handling
 Instrument Availability	Limited due to longer test time	Higher availability, fewer instruments needed

Benefits	Other Rapid Systems	ATLAS Rapid 20-Minutes System
 Test Duration	24 hours	20 minutes
 Result Receipt	Varies (manual/digital)	Automated Printout, minimizing human error
 Ease of Use	Varies	Touch screen display for simple operation and training
 Contamination Control	Varies	Lid-cover for spore tests to reduce contamination
 Cost	Standard	30% less expensive than leading brands
 Efficiency	High	Even higher efficiency, more emergency patient handling
 Instrument Availability	High	Higher due to quicker test time

COMBO! LABELLER & LABELS

ATLAS SHIELD LABELS



ATLAS SHIELD LABELLER

25% OFF



CLICK OR SCAN THIS QR Code to see Videos on Use

LABEL ROLLS
REG \$19.99/ROLL OF 1000
NET \$14.99/ROLL

REORDER CODE **PH-AT1000**

REG **\$299.98/UNIT**
PROMO
\$225.98/UNIT

REORDER CODE **PH-AT5600**

ATLAS SHIELD™ BIOLOGICAL INDICATORS - 24H

BOX OF 100

SAVE \$250!

NET \$166.63/BOX

2+1

REG \$249.95/BOX OF 100

REORDER CODE

AT-2492



2+1

BOWIE-DICK TEST CARDS

BOX OF 100

SAVE \$400!



NET \$174.99/PACK OF 100

2+2

REG \$349.99/PACK OF 100

REORDER CODE

AX-BD100



2+2

ATLAS SHIELD™ TYPE 5 CHEMICAL INTEGRATORS

PACK OF 250

5+5


NET \$19.99/PACK

5+5

REG \$39.95/PACK OF 250

REORDER CODE

AT-AXCP90



ATLAS WASH - U - TEST ULTRASONIC INDICATOR

PACK OF 50

2+2

NET \$74.97/PACK

2+2

REG \$149.95/PACK OF 50

REORDER CODE

PH-AT-50



\$49.99/HOLDER

FMD - STERILIZATION POUCHES CLASS 4 SELF-SEALING

200/BOX - DUAL INDICATOR

PROMO 6+4

#REORDER CODE	#SIZE	REG.	NET.
PH-P27510	2.75" X 10"	\$12.99	\$7.79
PH-P3510	3.5" X 10"	\$12.99	\$7.79
PH-P35625	3.5" X 6.25"	\$11.99	\$7.19
PH-P52511	5.25" X 11"	\$24.99	\$14.99

HB - STERILIZATION POUCHES (FINAL SALE)

CLASS 4
200 POUCHES/BOX

PROMO 10+10

#REORDER CODE	#SIZE	REG.	NET.
GSM442090	2.75" X 9"	\$9.99	\$4.99
GSM442101	3.5" X 5.25"	\$9.99	\$4.99
GSM442100	3.5" X 6.25"	\$11.99	\$5.99
GSM442102	5.25" X 11"	\$19.99	\$9.99
GSM441016	8" X 13"	\$29.99	\$14.99

FMD BARRIER FILM ROLL CLEAR & BLUE



NET \$17.99/ROLL

PROMO 3+1

REG. \$23.99/ROLL

CLEAR

BLUE

REORDER CODE

PH-BF02C

REORDER CODE

PH-BF01B

OPTIM OS1 WIPE 1 MINUTE KILL TIME 160/CAN



REORDER CODE

PH-OSICA-W12

INTRODUCING NEW HALYARD WRAPS SMS STERILIZATION WRAP



SIZE	BOX	PRICE	#REORDERS	10% DISCOUNT
12X12	500	\$79.95	PH-D31212	\$71.95
15X15	500	\$116.49	PH-D31515	\$104.84
20X20	500	\$168.95	PH-D32020	\$152.05
24X24	500	\$224.99	PH-D32424	\$202.49

INTRODUCING HB CSR WRAP



SIZE	BOX	PRICE	#REORDERS	10% DISCOUNT
12X12	500	\$74.99	PH-HBCSR1212	\$67.49
15X15	500	\$89.99	PH-HBCSR1515	\$80.99
20X20	500	\$149.99	PH-HBCSR2020	\$134.99
24X24	500	\$198.99	PH-HBCSR2424	\$179.09

OPTIM® 33 TB SURFACE WIPES

CAN OF 160 WIPES

BUY 24

BUY 12



REORDER CODE

PH-SCI33W12

BUY 24

\$16.95/CAN OF 160 WIPES

BUY 12

\$19.95/CAN OF 160 WIPES

OPTIM 33TB LIQUID 4L

BUY 4



REG **\$59.99**/BTL

BUY 4

NET **\$49.49**/BTL

REORDER CODE

PH-SCI334X4

MICRYLIUM BIO MERS

ASEPTIC TRANSPORT SOLUTION
5L BAG IN BOX

SAVE
10%



NET **\$89.93**/BX

10% DISCOUNT

REG **\$99.95**/1 BX

REORDER CODE

PH-02-MERS-005

INSTRUMENT CLEANER!

SCICAN

OPTIM CS (COLD SOAK),
CHEMOSTERILANT, 4L

BUY 4



NET **\$59.90**/BTL

BUY 4

REG **\$69.99**/1 BTL

REORDER CODE

PH-SCICS4X4

MICRYLIUM LINES - 5L

5L BAG IN BOX

SAVE
10%



NET **\$109.95**/BX

10% DISCOUNT

REG **\$119.95**/1BX

REORDER CODE

PH-03-BLUE-005

BIO TEXT 5L

CLINICAL SURFACE
DISINFECTANT

10%
DISCOUNT



NET **\$89.95**/BOX

10% DISCOUNT

REG **\$99.95**/1BX

REORDER CODE

PH-01-TEXT-005

STELLA® 300 NITRILE EXAM GLOVES

POWDER FREE | LATEX FREE | NON STERILE

XS S L

BOX OF 300

10+10



NET \$17.49/BOX

10+10

REG \$34.99/BOX OF 300

	XS	S	L
REORDER CODE	PH-ADV100SB	PH-ADV100MB	PH-ADV100LB
REORDER CODE	PH-ADV100LB	PH-ADV100MB	PH-ADV100XB

HALYARD AQUASOFT NITRILE GLOVE

XS S M L

BOX OF 300

7+3



NET \$22.39/BOX

7+3

REG \$31.98/BOX OF 300

	XS	S	M	L
REORDER CODE	PH-HA-300XS	PH-HA-300S	PH-HA-300M	PH-HA-300L

INTCO BLACK NITRILE GLOVES

S M L XL

BOX OF 100

6+4



NET \$7.71/BOX

6+4

REG \$12.85/BOX OF 100

	S	M	L	XL
REORDER CODE	PH-ADV100SB	PH-ADV100MB	PH-ADV100LB	PH-ADV100XB

RONCO SILKTEX LATEX EXAM GLOVES - 5MIL

XS S M L

BOX OF 100

6+4



NET \$10.19/BOX

6+4

REG \$16.99/BOX OF 100

	XS	S	M	L
REORDER CODE	PH-MT-819XS	PH-MT-819S	PH-MT-819M	PH-MT-819L

INTCO NITRILE GLOVES

S M L XL

BOX OF 100

10+10



NET \$5.25/BOX

10+10

REG \$10.50/BOX OF 100

	S	M	L	XL
REORDER CODE	SY-100S	SY-100M	SY-100L	SY-100XL

PRI-MED BLUEBERRY NITRILE GLOVES

S M L XL

BOX OF 250

7+3



NET \$20.99/BOX

7+3

REG \$29.99/BOX OF 200

	S	M	L	XL
REORDER CODE	PM6-2806	PM6-2807	PM6-2808	PM6-2809



ADELLE 3

LEVEL 3 ASTM MASKS

\$7.99 NET **7+3**
\$13.69

ALUMINUM NOSE CLIP
FOR STRONG HOLD

ULTRA SOFT
&
COMFORTABLE

● BLACK

REORDER CODE

PHB-180

50/BOX

MADE SPECIALLY FOR DENTAL PROFESSIONALS

ADELLE™ MASKS

Non-woven Fabric | Provides comfort & Protection | Single Use

ADELLE 3

ANTI-FOG

LEVEL 3 ASTM MASK

50/BOX

NET \$11.25/BOX

9+3

REG **\$14.99/BOX**

REORDER CODE

PHB-180A



BRONX AIR WATER SYRINGE TIPS

PACK OF 250

4+2



NET \$29.94/PACK

4+2

REG \$44.90 /PK OF 250

REORDER CODE

BRO-2381

ZEN SALIVA EJECTORS CLEAR WITH WHITE TIP

PACK OF 100

8+2



NET \$4.76/PACK

8+2

REG \$5.95 /PK OF 100

REORDER CODE

ZEN-2384

CURVED UTILITY #412 SYRINGE 12CC

PACK OF 50

BUY 4



NET \$29.99/BAG

BUY 4

REG \$34.99 /PK OF 50

REORDER CODE

100-5212

LEO DENTAL BIBS BLUE & PURPLE

PACK OF 500

3+2



NET \$26.99/CASE

3+2

REG \$44.99 /PK OF 250

BLUE

REORDER CODE

LEO-2341

PURPLE

REORDER CODE

PH-LEOLAV

REN™ TRAY COVER BLUE

CASE OF 1000

3+1



NET \$24.74/CASE

3+1

REG \$32.99 /CASE OF 1000

REORDER CODE

REN-2386

HB TIPS MIXING TIPS

BLUE - 1:1 - REORDER CODE - 100-1418
PINK - 5.0MM - 100-1417 - REORDER CODE - 100-1417
YELLOW - 4.2MM - REORDER CODE - M3100-1412
PURPLE - 7.5MM - REORDER CODE - M3100-1416
M3 HIGH PERFORMANCE MIXING TIPS
(GREEN - 6.5MM - HEAVY BODY & BITE REGISTRATION)
- REORDER CODE - PH-100-1413

NET \$29.69/BOX

10% DISCOUNT

REG \$32.99 /BOX

10%
DISCOUNT



CHEETA PROPHY CUPS

LATCH TYPE

CHE-PCLS-BL SOFT - BLUE
CHE-PCLS-PI SOFT - PINK

PACK OF 144

3+1



NET \$26.24/PACK

3+1

REG \$34.99 /PACK

REORDER CODE

PH00019

REORDER CODE

PH-CHE-PCLS-PI

SCREW TYPE

CHE-PCS-BL SOFT - BLUE
CHE-PCS-PI SOFT - PINK

PACK OF 144

3+1



NET \$24.74/PACK

3+1

REG \$32.99 /PACK

REORDER CODE

PH-CHE-PCS-BL

REORDER CODE

PH00020

SNAP-ON

CHE-PCSOS-B SOFT - BLUE

PACK OF 144

3+1



NET \$14.99/PACK

3+1

REG \$19.99 /PACK

REORDER CODE

PH-CHE-PCSOS-B

MILO

PRE-BENT APPLICATOR TIPS

18 GAUGE - PINK | REORDER | MIL-2387
20 GAUGE - BLACK | REORDER | MIL-2388
22 GAUGE - GREY | REORDER | MIL-2389
25 GAUGE - BLUE | REORDER | MIL-2390

NET \$11.19/BAG

4+1

REG \$13.99 /BAG OF 100



CHEETA HEADREST COVER

CASE OF 500 - 10X13

3+1



NET \$40.49/CASE

3+1

REG \$53.99 /CASE

BLUE

REORDER CODE

PH00011

WHITE

REORDER CODE

PH00010

KAI MICRO APPLICATORS

REGULAR | REORDER | KAI-2304
FINE | REORDER | KAI-2302
SUPERFINE | REORDER | KAI-2301

3+1



BOX OF 400

NET \$26.21/BOX

3+1

REG \$34.95 /BOX OF 400
4 TUBES OF 100

MARK 3

INTRAORAL MIXING TIPS

INTRAORAL TIPS YELLOW | REORDER | 100-0060
INTRAORAL TIPS CLEAR | REORDER | 100-0059

20%
DISCOUNT

PACK OF 100

NET \$27.99/PACK

20% DISCOUNT

REG \$34.99 /PACK OF 100



MARK3 GENERAL PURPOSE ULTRASONIC SOLUTION 4L



NET \$35.99/JUG

20% DISCOUNT

REG \$44.99/JUG

REORDER CODE

PH-100-0092

BIOENZYMATIC ULTRASONIC CLEANER 83036 64 TABLETS/BOX



NET \$51.99/BOX

4+1

REG \$64.99/BOX

REORDER CODE

M3100-7640

BM6000 EVACUATION SYSTEM CLEANER 4L



NET \$58.49/JUG

10% DISCOUNT

REG \$64.99/JUG

REORDER CODE

PH01680

HELIX™ T-STYLE LIGHT HANDLE COVER

BOX OF 500



NET \$14.96/BOX

3+1

REG \$19.95/BOX OF 500

REORDER CODE

HEL-2308

CHEETA TRAY SLEEVES PLASTIC 10.5"X14"

BOX OF 500



NET \$29.99/BOX

BUY 3+

REG \$36.99/BOX

REORDER CODE

PH00035

HB - HALF CHAIR COVER ON ROLL

27-1/2" X 24" 700MM X 610MM

BOX OF 225



NET \$44.99/BOX

10% DISCOUNT

REG \$49.99/BOX

REORDER CODE

HA-5011

HB - FULL CHAIR COVER ON ROLL



NET \$58.49/ROLL

BUY 3+

REG \$64.99/ROLL

REORDER CODE

FS-5012

MONOJECT DENTAL NEEDLES - PLASTIC

BOX OF 100



25 GAUGE LONG - RED 1-1/4" | REORDER CODE | MDT8881400033
25 GAUGE SHORT - RED 3/4" | REORDER CODE | MDT8881400041
27 GAUGE LONG - PLASTIC HUB - YELLOW 1-1/4" | REORDER CODE | MDT8881400058
30 GAUGE EXTRA SHORT 1/2" BLUE | REORDER CODE | MDT8881400073
27G SHORT 3/4" - PLASTIC HUB | REORDER CODE | MDT8881400056
30 GAUGE SHORT - PLASTIC HUB - BLUE 3/4" | REORDER CODE | MDT8881400074

NET \$18.74/BOX

3+1

REG \$24.99/BOX OF 100 NDL

CHEETA

NON WOVEN SPONGE - GAUZE - 2X2

CASE OF 5000

4+2



NET \$39.96/CASE

4+2

REORDER CODE

PH0007

REG \$59.95/CS OF 5000

ODIN

NON WOVEN SPONGE - GAUZE - 4X4

CASE OF 2000

4+2



NET \$43.32/CASE

4+2

REORDER CODE

ODI-2352

REG \$64.99/CS OF 2000

HUNTER™

PREMIUM DENTAL ROLL - SIZE 2

CASE OF 2000

2+1



NET \$23.30/CASE

2+1

REORDER CODE

HUN-2361

REG \$34.95/CS OF 2000

CHEETA

PAPER CUP CASE OF 1000

BUY 3+



NET \$39.99/CASE

2+1

REORDER CODE

PH-CHE-PC50-PC

REG \$46.99/CS OF 1000

AVA™

COTTON TIPPED APPLICATORS
BOX OF 1000



6"

3"

NET \$7.46/CASE

3+1

REORDER CODE - 3"

AVA-2313

REORDER CODE - 6"

AVA-2312

REG \$9.95/CS OF 1000

DISPOSABLE

EVACUATION TRAPS

STARTING FROM \$64.95

SIZES & REORDER NUMBERS

2-3/4" X 3-5/8" | REORDER | 100-12200 | \$94.99

3-1/2" X 4-3/8" | REORDER | 100-12300 | \$94.99

1-1/4" | REORDER | 100-5500 | \$64.95

2-1/8" - 5501 | REORDER | 100-5501 | \$64.95

1-7/8" | REORDER | 100-5503 | \$64.95

3-3/4" | REORDER | 100-5505 | \$64.95

1.5" | REORDER | 100-5509 | \$64.95



4D RUBBER

DAM ISODAM NON-LATEX

SIZES & REORDER NUMBERS

5" X 5" - MEDIUM BLUE | REORDER | 4DM-800520

5" X 5" - HEAVY BLUE | REORDER | 4DM-300520

6" X 6" - MEDIUM BLUE | REORDER | 4DM-800615

6" X 6" - HEAVY BLUE | REORDER | 4DM-300615



NET \$31.19/BOX

4+1

REG \$44.99/BOX

4+1

MONOJECT DENTAL NEEDLE - METAL

COVIDIEN MONOJECT DENTAL NEEDLES - 100

400 METAL HUB

25G GAUGE RED 1-1/4" SHORT | REORDER | PH01820

25G GAUGE RED LONG | REORDER | PH01818

27G GAUGE YELLOW 1-1/4" - SHORT | REORDER | MDT8881401064

27G GAUGE YELLOW LONG - SHORT | REORDER | PH-Z8881401056

30G GAUGE BLUE SHORT | REORDER | PH-Z8881401072



REG \$26.99/BOX

NET \$21.99/BOX



LOOK SUTURES

REVERSE CUTTING NEEDLE - 12/BOX

\$48.99/BOX

PLAIN GUT REVERSE CUTTING

551B C6 4-0 27"	L00551B
553B C7 3-0 27"	L00553B
554B C6 3-0 27"	L00554B
555B C7 2-0 27"	L00555B
591B C6 4-0 18"	L00591B
592B C6 3-0 18"	L00592B
594B C31 3-0 18"	L00594B
596B C3 5-0 18"	L00596B

CHROMIC GUT REVERSE CUTTING

525B C6 3-0 18"	L00525B
544B C26 4-0 18"	L00544B
557B C6 5-0 27"	L00557B
558B C6 4-0 18"	L00558B
559B C6 4-0 27"	L00559B
560B C6 3-0 27"	L00560B
563B C17 5-0 18"	L00563B
595B C3 4-0 18"	L00595B

TOPICAL ANESTHETIC GEL

BUBBLE GUM GEL 10Z | REORDER | M3100-1580
 MINT GEL 10Z - 1584 | REORDER | M3100-1584
 PINA COLADA GEL 10Z | REORDER | M3100-1585
 CHERRY GEL 10Z | REORDER | M3100-1581
 STRAWBERRY GEL 10Z | REORDER | M3100-1586



NET \$9.90/EACH

BUY 3

REG \$12.95/EACH

ENDO ICE



NET \$35.99/CAN

4+1

REG \$44.99/CAN

REORDER CODE

PH01836

KEYSTONE PROPHY PASTE - 340G

1.23% FLUORIDE ION GELATO

JAR - 340ML/12OZ - MINT - COARSE | REORDER | PH-24-04077
 340 ML MEDIUM MINT | REORDER | KEY24-08177M
 340ML/12OZ - BUBBLE GUM - FINE | REORDER | PH-24-10877



NET \$21.90/EACH

BUY 3

REG \$24.99/BOX

VPS IMPRESSION MATERIAL

Light Body, Medium Body Available in Reg or Fast Set
 4 X 50ML CARTRIDGES

LIGHT BODY - REGULAR SET - 4 X 50ML CARTRIDGES | REORDER | 100-3010
 LIGHT BODY - FAST SET - 4 X 50ML CARTRIDGES | REORDER | 100-3011
 MEDIUM BODY - REGULAR SET - 4 X 50ML CARTRIDGES | REORDER | 100-3012
 MEDIUM BODY - FAST SET - 4 X 50ML CARTRIDGES | REORDER | 100-3013
 HEAVY - REGULAR SET - 4 X 50ML CARTRIDGES | REORDER | PH01985
 HEAVY BODY - FAST SET - 4 X 50ML CARTRIDGES | REORDER | 100-3015

10% Off

REG \$59.99/BOX

NET \$53.99/BOX

BITE REGISTRATION MATERIAL

2 X 50ML CART. + MIXING TIPS
 REGULAR SET | REORDER | PH-0218
 FAST SET | REORDER | 100-0219

10% Off

REG \$49.99/BOX

NET \$44.99/BOX

VPS PUTTY

300 ML EACH BASE AND CATALYST
 REGULAR SET 300ML BASE & 300ML CATALYST | REORDER | 100-3010
 FAST SET | REORDER | 100-3011

10% Off

REG \$79.99/BOX

NET \$71.99/BOX

VPS TRAY ADHESIVE

10 ML BOTTLE
 REORDER | PH-100-3000



NET \$20.69/BOX

10% Off

REG \$22.99/BOX

MARK 3 APF FLUORIDE GEL 17OZ

17OZ MINT | REORDER | PH-100-1612
 17OZ BUBBLE GUM | REORDER | PH-100-1613
 STRAWBERRY 17OZ. | REORDER | 100-1614



NET \$19.99/EACH

BUY 3+

REG \$24.99/EACH



ALGINATE ALTERNATIVE

(VINYL POLYSILOXANE)

FAST SET - 50ML CARTRIDGE - 8 CARTRIDGES - MINT FLAVOUR - PINK COLOUR - COMPARES TO ALGINOT

WORK TIME: 1 MINUTE
 INTRA ORAL SET TIME: 1 MINUTE



NET \$119.90/EACH

BUY 2

REG \$149.95/EACH

REORDER CODE

PH-100-9010

DIA ROOT

BIO SEALER

- BioCeramic: excellent biocompatibility
- pH 12: highly antibacterial
- Superior bonding to dentin & GP
- Low solubility: zero shrinkage
- Washout resistant
- Optimal flowability of 22mm
- Saves time: pre-mixed and pre-filled

#1003-701 Dia-Root Bio Sealer
 Contains: 1 x 2g syringe, 20 tips, 1 cap

NET \$96.74/EACH

BUY 3+1

REG \$128.99/EACH

REORDER CODE

PH-1003-701



DIA ETCH

#A2001-3101 ECONO PACKAGE

CONTAINS: 5 SYRINGES X 5 ML + 20 DISPOSABLE TIPS



NET \$29.99/EACH

BUY 3+2

REG \$49.99/EACH

REORDER CODE

A2001-3103

PT-A™ | AIR POLISHER & ULTRASONIC SCALER

Elevate your practice today!



CLICK OR SCAN THIS QR CODE FOR VIDEOS AND RESOURCES!

0%

Financing For 24 months

\$332.98/MONTH

*terms and conditions apply

PROMO
GET 25 POWDER BOTTLES FREE

REORDER CODE

PH01992




PROPHYLAXIS POWDERS

FAST

Sodium Bicarbonate
(Supra - 65um)

GENTLE MINI

Glycine
(Sub & Supra - 25um)



Recommended angle:
30°-60°



Recommended distance:
3-5mm

MOBILE
CART

SPECIAL
\$798/UNIT

REORDER CODE

PH-WP-MC



WOODPECKER PT-A / PT-B

WOODPECKER PT-A / PT-B
SUPRAGINGIVAL HANDPIECE
(AP-1) - AP-1

3+1

NET \$449.96/UNIT
(MIX AND MATCH)

3+1

REG **\$599.95**/UNIT



REORDER CODE

PH-WP-AP-1

WOODPECKER PT-A / PT-B
SUBGINGIVAL HANDPIECE
(AP-2) - AP-2

3+1

NET \$449.96/UNIT
(MIX AND MATCH)

3+1

REG **\$599.95**/UNIT



REORDER CODE

PH-WP-AP-2

WOODPECKER PT-A / PT-B
ULTRASONIC HANDPIECE
(TZ-2L) - TZ-2L

NET \$449.96/UNIT
(MIX AND MATCH)

3+1

REG **\$599.95**/UNIT



REORDER CODE

PH-WP-TZ-2L

WOODPECKER PT-A / PT-B
SUBGINGIVAL AIR POLISHING
NOZZLES

2+1

NET \$46.63/PACK

2+1

REG **69.95**/PK OF 10 NOZZLES



REORDER CODE

PH-WPNOZZLE10

WOODPECKER DTE | ULTRASONIC SCALING TIPS

NET \$69.95/TIP

P33 TIP



REORDER CODE

PH-P33

P50L TIP



REORDER CODE

PH-P50R

P50R TIP



REORDER CODE

PH-P59

P56 TIP



REORDER CODE

PH-P50L

P59 TIP



REORDER CODE

PH-P56

E60 TIP



REORDER CODE

PH-E60

2+2

NET \$69.95/TIP

REG **\$139.95**/TIP

WOODPECKER GENTLE

(GLYCINE - SUPRAGINGIVAL 65UM - 200G BOTTLE - PROPHYLAXIS POWDER)

NET \$24.99/BTL

PH-WP-GENTLE

2+2

REG **\$49.99**/120G BTL



WOODPECKER GENTLE MINI

(GLYCINE - SUBGINGIVAL + SUPRAGINGIVAL POWDER 25UM - 120G BOTTLE - PROPHYLAXIS POWDER)

NET \$24.99/BTL

PH-WP-GENTLE-MINI

2+2

REG **\$49.99**/120G BTL



WOODPECKER FAST

(GLYCINE - SUPRAGINGIVAL 65UM - 200G BOTTLE - PROPHYLAXIS POWDER)

NET \$24.99/BTL

PH-WP-FAST

2+2

REG **\$49.99**/120G BTL





LUNA
WHITENING

SEND YOUR PATIENTS HOME WITH A GLOWING SMILE!

LUNA WHITENING TAKE HOME KIT

6+4

A NEW ADVANCED TAKE HOME TOOTH WHITENING SYSTEM WITH LED LIGHT TECHNOLOGY

INCLUDED

- 2 WHITENING SYRINGES X 5ML
- ACTIVATION LIGHT
- SILICONE TRAY
- GUIDE & INSTRUCTIONS



NET \$54.97

5+5

REG **\$109.95** /KIT

REORDER CODE

PH-LUNATHK

LUNA/PURE

GINGIVAL (GUM) BARRIER SYRINGE
(3ML - TREATS 3 CLIENTS)

2+1

NET \$33.32 /SYRINGE

2+1

REG **\$49.99** /SYRINGE



REORDER CODE

PH-LGBS3

LUNA WHITENING SMILE BOOST PACK

2 WHITENING SYRINGES X 3ML - 35% CP

4+2

NET \$19.99 /KIT

4+2

REG **\$29.99** /PACK

REORDER CODE

PH-LTH2



LUNA DESENSITIZING GEL (5ML SYRINGE)

- STRENGTHENS AND REMINERALIZES TOOTH ENAMEL, REDUCING SENSITIVITY AFTER WHITENING.
- PROMOTES ACTIVE REMINERALIZATION WITH FLUORIDES, CALCIUM, AND PHOSPHATE.
- EASY TO APPLY, EFFECTIVE FOR NIGHTTIME USE.
- CONTAINS SODIUM FLUORIDE TO BIND MINERALS TO ENAMEL.

1+1

NET \$14.99 /SYRINGE

1+1

REG **\$29.99** /SYRINGE

REORDER CODE

PH-DG5S



PROFESSIONAL ADVICE

Dental Professionals - Select the right LUNA Whitening Concentration by assessing the patient's oral health and sensitivity. Use a lower concentration (16% CP to 25% CP) for sensitive teeth and higher concentrations (35% HP) for faster results. Tailored treatments ensure safe, effective, and satisfiable whitening! (Devila et al, 2020)*

ACERO STAINLESS STEEL CROWNS

ACERO NARROW

Acero Narrow Molar Crowns



Less mesial distal width to treat that large class 2 caries with space loss.

ACERO
CROWNS
www.acerocrowns.com

NET \$36.66/BX

20+10

REG \$54.99/BX



ACERO CROWN ADVANTAGE


- ⊙ Simplified universal numbering system crown label correspond to tooth number
- ⊙ Surgical-grade 316 stainless steel
- ⊙ User Friendly tray layout
- ⊙ Secure traycase

ACERO CROWNS^{3S}

Acero3S Posterior Crowns



- Laser-etched sizes (autoclave-proof = less waste!)
- Made in the U.S. by veterans
- 100% compatible and interchangeable with the leading brands
- Preformed and pre-contoured to prevent microgaps
- Acero 3S molar crowns available in the standard 48 sizes
- tested by a leading US dental school and shown to have superior fracture resistance.



ACERO CROWNS^{XT}

AceroXT Original Crowns

- Acero XT original all metal crowns are available in the standard 48 molar sizes and 36 anterior sizes.

Developed by Dentists,
Inspired by Years
of Clinical Experience



Trusted, esthetic and affordable dental crowns endorsed by top dentists for making your kid's perfect smile last forever.

NET \$36.66/BX

20+10

REG \$54.99/BX

Why Acero Crowns?

Acero Crowns are crafted from surgical-grade 316L stainless steel more resistant to corrosion and wear than 304 stainless steel used by competitors. The crowns are preformed and pre-contoured (or commonly referred to as "pre-cripped") with a shallow occlusal anatomy to preserve healthy tooth structure with less tooth reduction.

Competitor Comparison Chart

FEATURE	ACERO CROWNS 3S	COMPETITION	DESCRIPTION
STEEL GRADE	*316L Grade	304 Grade	Acero3S uses surgical-grade stainless steel which is more resistant to corrosion/wear
NUMBERING	*Laser Etched	Ink	Acero3S sizes are laser-etched and won't wash away after being autoclaved
SAVINGS	*40%-65%	0%	Acero3S crowns cost 40%-65% less than other leading brands' stainless steel crowns
SHAPE	Bell	Bell	This shape accurately duplicates the anatomy of primary molars
CRIMPED	Cervical Margin	Cervical Margin	Provides for proper retention and a snap fit
PREFORMED	Life-like Height	Life-like Height	Allows for appropriate contour and occlusal surface
TRIMMING	Yes	Yes	Trimming can be performed with crown scissors or an abrasive wheel
ADHESIVE	Multiple	Multiple	Resin-modified glass ionomer, polycarboxylate or zinc phosphate

Acero Crowns



Half the price of other leading brands

Laser-etched sizes (autoclave-proof=less waste!)

Surgical-grade 316L U.S stainless steel

100% compatible and interchangeable with the leading brands

Preformed and pre-contoured to prevent microgaps

**ACERO
CROWNS**



PRECISION EDGE SYSTEMS SPECIALS

COMPARES TO PALODENT SYSTEM

PENGUIN PRECISION EDGE KIT SECTIONAL MATRIX SYSTEM KIT

CONTAINS:

100 MATRICES - 25 EACH
3.5 mm - 4.5 mm - 5.5 mm - 6.5 mm

75 WEDGES - 25 EACH

Small - Medium - Large

- 1 UNIVERSAL RING
- 1 NARROW RING
- 1 FORCEPS & 1 PIN TWEEZERS



NET \$359.40

INTRO KIT A

40% OFF

REG **\$599**/KIT

REORDER CODE

PH-EDGEKIT

PENGUIN PRECISION EDGE RING REFILL

COMPARES TO PALODENT RING
MADE WITH NITI MATERIALS

- 1 UNIVERSAL RING | PH-DMB0495
- 1 NARROW RING | PH-DMB0496

**MORE FLEXIBLE AND
DURABLE THAN
COMPETITOR!**



1 PER BOX

NET \$83.32/1BX

2 + 1

REG **\$124.99**/BX

PENGUIN PRECISION TULWAR MATRICE WEDGE REFILL

- CAN BE DOUBLE-WEDGED
- THE TIP ANGLED TO PROTECT RUBBER DAM
- EASY TO INSERT AND REMOVE

SMALL/DARK BLUE | PH-DIW0903
MEDIUM/LIGHT BLUE | PH-DIW0902
LARGE/WHITE | PH-DIW0901



Small



Medium



Large

**COMPARES TO PALODENT
MATRICE WEDGE**

BOX OF 100

NET \$44.99/BX OF 100

3 + 1

REG **\$59.99**/BX

PENGUIN PRECISION MATRICE REFILL

- SUBGINGIVAL MATRICES AVAILABLE!

**COMPARES TO PALODENT
MATRICE**

3.5MM | PH-DMB0811
PRE-MOLAR - 4.5MM | PH-DMB0812
MOLAR - 5.5MM | PH-DMB0813
MOLAR W/EXTENSIONS - 6.5MM | PH-DMB0814



3.5mm



4.5mm



5.5mm



6.5mm

BOX OF 50

NET \$52.49/BX OF 50

3 + 1

REG **\$69.99**/BX

PENGUIN PRECISION PRO SYSTEMS SPECIALS

COMPARES TO GARRISON SYSTEM

PENGUIN PRECISION PRO KIT SECTIONAL MATRIX SYSTEM KIT

CONTAINS:

100 MATRICES - 25 EACH

3.5 mm - 4.5 mm - 5.5 mm - 6.5 mm

75 WEDGES - 25 EACH

Small - Medium - Large

- 3** PRO RING (SMALL, MEDIUM, LARGE)
- 1** PIN TWEEZERS
- 1** FORCEP



NET \$539.40

INTRO KIT A

40% OFF

REG **\$899** /KIT

REORDER CODE

PH-PRECISIONPROKIT

PENGUIN PRECISION PRO RING REFILL

COMPARES TO GARRISON RING
MADE WITH NITI MATERIALS

SMALL/BLUE | PH-DMB0490

MEDIUM/ORANGE | PH-DMB0491

LARGE & WIDE/GREEN | PH-DMB0494

**MORE FLEXIBLE AND
DURABLE THAN
COMPETITOR!**

1 PER BOX



NET \$86.66 /BX OF 1

2+1

REG **\$129.99** /BX

PENGUIN PRECISION WEDGE GUARDS

MATRICE WEDGE REFILL

- EASY TO INSERT AND REMOVE

**COMPARES TO PALODENT
WEDGE GUARDS**

SMALL - YELLOW | PH-DMB0508

MEDIUM - ORANGE | PH-DMB0507

LARGE - PURPLE | PH-DMB0506



BOX OF 48



NET \$52.49 /BX OF 48

3+1

REG **\$69.99** /BX

PRECISION PRO EZY CONTACT

MATRICE REFILL

- NON-STICK COATING

**COMPARES TO PALODENT &
GARRISON**

3.5MM | PH-DMB1801

PRE-MOLAR - 4.5MM | PH-DMB1802

MOLAR - 5.5MM | PH-DMB1803

MOLAR W/EXTENSIONS - 6.5MM | PH-DMB1804



BOX OF 50



NET \$67.49 /BX OF 50

3+1

REG **\$89.99** /BX

COMPARES TO SOFLEX DISCS

PRECISION

POP-ON COMPOSITE FINISHING DISCS

- COLOUR CODED
- ADAPTIVE
- EASY MANDREL ATTACHMENT



3+1

REFILL 1/2"
85/PACK

Course	Medium	Fine	Super Fine
(M1053)	(M1054)	(M1055)	(M1056)
PH-M1053	PH-M1054	PH-M1055	PH-M1056

REFILL 3/8"
85/PACK

Course	Medium	Fine	Super Fine
(M1057)	(M1058)	(M1059)	(M1060)
PH-M1057	PH-M1058	PH-M1059	PH-M1060

NET \$52.49

3+1

REG \$69.99 /PK

NET \$46.66

6+3

REG \$69.99 /PK

MANDREL RA REFILL

3/ PACK
(PH01741)



NET \$35.99

4+1

REG \$44.99 /PK

PRECISION STEM POLISHING DISCS



MANDREL



STANDARD SHORT

REORDER CODE

PH-DPW0400

REORDER CODE

PH-DPW0400B

FULL KIT

100 DISCS (25/GRIT) +
4 MANDREL

\$69.99

INTRO KIT

40 DISCS (10/GRIT) +
2 MANDREL

\$39.99



STEM POLISHING DISCS



COARSE

MEDIUM

FINE

SUPER FINE

O-STAR 7 Modes

Unprecedented
Curing Light
incorporating
various strengths

Charging base with
built-in light meter

- 7 MODES (NORMAL, HIGH, TURBO, ORTHO, SOFT, PULSE, AND CARIES CHECKER)
- INDUSTRY-LEADING PERMEABILITY AND LIGHT-LOSS REDUCTION!
- SLIM-HEAD AND 360-DEGREE ROTATION
- WIDE-RANGE WAVELENGTH CURING!

BUY 20

MULTICHROME COMPLUES BOXES

GET O STAR CURING LIGHT FREE
SAVINGS OF

\$1300

REORDER CODE

PH-WP-O-STAR





HARVARD®

Harvard **MultiChrome**

NEW



**It's
magic.**



EN

Harvard MultiChrome

It's magic.

Shade-adaptive
composite for almost all
tooth shades.

- Hyper-Nano technology, spherical fillers
- Without any color pigments
- Light-optical chameleon effect for perfect match with natural tooth substance
- Simple storage logistics



The Next



CLICK OR SCAN QR CODE FOR
BROCHURES & RESEARCH!



MULTICHROME UNIVERSAL COMPULES

20 COMPULES X
0.25G/PACK!



NET \$134.99 /PK

12+4

REORDER CODE

PH01855

NET \$143.99 /PK

4+1

REG \$179.99 /PACK

MULTICHROME UNIVERSAL 3G SYRINGE



NET \$112.49 /PK

12+4

REORDER CODE

PH01854

NET \$119.99 /PK

4+1

REG \$149.99 /PACK

MULTICHROME BLOCK N MASK COMPULES

BUY MULTICHROME
AND GET SPECIAL PRICE!

BOX OF 20 COMPULES



NET \$134.99 /PK

12+4

REORDER CODE

PH-7100007

NET \$143.99 /PK

4+1

REG \$179.99 /PACK

MULTICHROME UNIVERSAL FLOW SYRINGE PACK



NET \$97.46 /PK

12+4

REORDER CODE

PH-7100003

NET \$103.96 /PK




4+1

REG \$129.95 /PACK

Harvard UltraFill

Ultra Gloss, Opalescent, High Performance

Composite for anterior and posterior restorations.



Harvard UltraFill

Ultra Gloss, Opalescent, High Performance Composite
for **anterior** and posterior restorations

- Permanent high gloss
- Opalescence
- Tooth-like fluorescence
- Easy to work with
- Superior handling properties
- Low polymerization shrinkage
- High abrasion resistance
- Optimal physical properties
- Radiopaque

seit **1892**
HARVAR



IntroKit S, 3 x 3 g



3 g



20 x 0.25 g

HARVARD ULTRAFILL COMPULES

NET \$142.49 /PK	12+4	
NET \$151.96 /PK	4+1	REG. \$189.99 /PACK

HARVARD ULTRAFILL SYRINGE SPECIALS

NET \$112.46 /PK	12+4	
NET \$119.96 /PK	4+1	REG. \$149.95 /PACK

Universal Opacity easy technique



U

Universal



CLICK OR SCAN QR CODE FOR
BROCHURES & RESEARCH!

- Selected shades in **one** opacity
- In 80 % of clinical cases best compromise of opacity

Harvard UltraFill „Universal Opacity“		
Shade	3 g syringe	20 x 0.25 g OptiTips®
A1 U	7111111	7111211
A2 U	7111112	7111212
A3 U	7111113	7111213
A3.5 U	7111114	7111214
BL U	7111151	
Harvard Applier OptiTips®		7095200

Harvard InterLock® ONE Universal Adhesive

Really universal: suitable for all tooth conditioning techniques:

as "Self Etch" adhesive, as "Total Etch" adhesive (which etching step), wet or dry

- Only one thin layer necessary
- High & reproducible bond strength to enamel and dentine even under long term stress
- Effective marginal seal for low risk of postoperative sensitivity
- Excellent performance also with self or dual cure composites when light curing is possible, if not see e.g. for adhesive post cementing in the root canal mix additional with Harvard InterLock ONE Self Cure Activator
- Convenient application by bottle (precise dropper)



Harvard InterLock® ONE – Universal Performance

Dentin & enamel **after thermocycling**

Etching: Total Etch or self etch

Composite: Light cured – **InterLock ONE**: Light cured

Composite: Self cured – **InterLock ONE**: Light cured

Composite: Self cured – **InterLock ONE** plus **InterLock ONE Self Cure Activator** ★ if light cure is not possible eg. in the root canal.




Article information	Order no.
Harvard InterLock ONE	
Universal Adhesive	
5 ml bottle	7083611

**HARVARD
INTERLOCK ONE
(5ML)**



REORDER CODE
7083611



NET **\$179.96/BTL**
BUY 4+1
REG **\$224.95/BTL**

CLICK OR SCAN QR CODE FOR
BROCHURES & RESEARCH!




Harvard ZirconCore

Dual cure core build-up
and post cementation composite



- Excellent consistencies for both indications **core build-up** and **post cementation**
- “Cuts like dentin“
- Very good mechanical properties for permanent restorations
- Contains a substantial amount of Nano-zirconia
- Fluoride release
- Thixotropic properties
- Very good radiopacity
- Self-mixing material with the Minimix syringe
- For optimal results, use with **Harvard InterLock® ONE / Harvard InterLock® ONE Self Cure Activator** or with **Harvard Bond SE Dual**



Excellent mechanical properties in all curing modes.

2+1

NET \$126.66

2+1



REG \$189.99

REORDER CODE - A2
PH-7083599


REORDER CODE - A3
PH-7083600



Physical Properties



ZC Harvard ZirconCore
RE Voco Rebuilda
FC Dentsply FluoroCore 2+
MC Ivoclar Vivadent MultiCore Flow
LC DMG LuxaCore Z Dual Automix



Post cementation and core build-up



1. Prepared and bonded root canal
2. Application into the root canal
3. Placement of the post
4. Light cured core build-up
5. Shaped core build-up

Technical data

Working time (23 °C / 73 °F)	1:30 min
Setting time (self cure)	3:30 min
Light cure (optional)	40 sec


ORDER INFORMATION

SHADE A2 - 5ML PACK

SHADE A3 - 5ML PACK



Liability is excluded for all printing errors and omissions. Before using our Harvard products, the respective directions for use should be noticed in every case. All measurements are internal measurements of Harvard Dental International.



NET \$799
REG \$1199

AIRPEAK™

MIGHTY POWER



ORDER INFORMATION



TITANIUM



GRIPMAX



QUASPRAY



SEAL HEAD SYSTEM



A1005
AIRPEAK™ X600-M/K, KAVO BACKEND

NET \$799/BOX
2+1
REG \$1199.99/BOX

REORDER CODE
PH-A1005



A1006
AIRPEAK™ X600-M/N, NSK BACKEND

NET \$799/BOX
2+1
REG \$1199.99/BOX

REORDER CODE
PH-A1006

AIRPEAK™

CLASSIC LOW SPEED HANDPIECE



2+1 REG. \$299.95/UNIT
NET \$199.97/UNIT

(A1010) G100-H1
• Latch Type

REORDER CODE
PH-G100-H1



2+1 REG. \$299.95/UNIT
NET \$199.97/UNIT

(A1011) G100-H2
• Push Button

REORDER CODE
PH-G100-H2

1:1 Direct Drive



2+1 REG. \$399.95/UNIT
NET \$266.63/UNIT

(A1008) G100-S1H1
• For CA burs (ø2.35) • Latch type
• Max Speed: 30,000 min

REORDER CODE
PH-A1008


1:1 Direct Drive



2+1 REG. \$399.95/UNIT
NET \$266.63/UNIT

(A1009) G100-S1H2
• For CA burs (ø2.35) • Clean Head System
• Push Button Chuck • Max Speed: 40,000 min

REORDER CODE
G100-S1H2



2+1 REG. \$399.95/UNIT
NET \$266.63/UNIT

(A1012) G100-ST2
• For HP burs (ø2.35) • Stainless Steel Body
• Clean Head System • Max Speed: 40,000 min

REORDER CODE
PH-G100-ST2

PROPHY HANDPIECE

4:1 Reduction



2+1 REG. \$499.95/UNIT
NET \$333.95/BOX

(A1017) G200-SC
• For screw-in cups and brushes • Stainless Steel Body
• 360° Rotation • w/o Spray • Paste Defense System
• Max Speed: 5,000 min

REORDER CODE
PH-G200-SC

4:1 Reduction



2+1 REG. \$499.95/UNIT
NET \$333.95/BOX

G200-LA
• For CA burs (ø2.35) • Push Button Chuck
• Stainless Steel Body • Clean Head System
• Ball Bearing • Max Speed: 5,000 min

REORDER CODE
PH-G200-LA



STERILIZED

7 DIAMOND

SINGLE-PATIENT USE DIAMOND

Single Patient Use Diamond burs, are individually wrapped, Gamma-Sterilized, 25 burs per package. Manufactured with natural diamonds which provide more cutting surfaces than synthetic diamonds; also Swiss multi-coating technique is used to provide cutting efficiency with durability and less clogging. There are around 200 shapes and grits, one low price for all these burs.

NET \$78.99

BUY 2+

REG \$84.99

BOX OF 25

ROUND END TAPER

ORDER/ISO #	1-850.010	1-856.012	1-850.012	1-S855.014	1-856.014	1-850.014	1-S855.016	1-856.016	1-850.016	1-856.018
GRITS	C M F	C M	C M	SC	C M F	SC M F	SC M	SC M F XF	SC M F	SC F
Head/Lgh/Tip(mm)	1.0/10/0.6	1.2/7/0.7	1.2/10/0.7	1.4/6/0.9	1.4/8/0.9	1.4/10/0.9	1.6/6/1.1	1.6/8/1.1	1.6/10/1.1	1.8/8/1.4
Neo. #		1112.7	1112.10	1114.6	1114.8	1114.10	1116.6	1116.8	1116.10	1118.7

ROUND END TAPER

ORDER/ISO #	1-856N.018	1-856L.018	1-856.021	1-S856.024	1-856.024	1-848.010	1-847.012	1-847L.014	1-847.014	1-847.016
GRITS	C	SC M F	SC	SC	SC F	C M	C	C	SC M F	SC M F
Head/Lgh/Tip(mm)	1.8/8/1.2	1.8/9/1.4	2.1/8/1.6	2.4/7/1.8	2.4/7/1.8	1.0/10/0.7	1.2/8/0.8	1.4/9/0.7	1.4/8/1.0	1.6/8/1.2
Neo. #	1118.8	1118.9	1121.8	2424CS	2424			0914.9	0914.8	0916.8

FLAT END TAPER

FLAT END TAPER

MODIFIED FLAT END TAPER

MODIFIED CHAMFER












ORDER/ISO #	1-848.016	1-847.018	1-848.018	1-S846KR.016	1-847KR.016	1-845KR.018	1-847KR.018	1-845KR.025	1-878K.012	1-879K.012	1-878K.014
GRITS	SC M	SC M	SC M	SC M	SC M F	SC M	SC	M	C M	C M	SC M F
Head/Lgh/Tip(mm)	1.6/10/1.1	1.8/8/1.4	1.8/10/1.4	1.6/6/1.1	1.6/8/1.1	1.8/4/1.4	1.8/8/1.4	2.5/4/2.0	1.2/8/0.8	1.2/10/0.8	1.4/8/1.0
Neo. #	0916.10	0918.8	0918.10	0816.6S	0816.8	0818.4	0818.8	2525	1712.8	1712.10	1714.8












MODIFIED CHAMFER












MOD. FLAT END CYL.











FLAT END CYLINDER













ORDER/ISO #	1-879K.014	1-S877K.016	1-878K.016	1-878K.018	1-879K.018	1-878K.021	1-837KR.014	1-S835.010	1-835.010	1-835.012
GRITS	C M	C M	SC M F	SC F	SC M	C	C F	C M	C M	C
Head/Lgh/Tip(mm)	1.4/10/1.1	1.6/6/1.2	1.6/8/1.2	1.8/8/1.4	1.8/10/1.4	2.1/8/1.6	1.4/8/1.4	1.0/4/1.0	1.0/4/1.0	1.2/4/1.2
Neo. #	1714.10	1716.6S	1716.8	1718.8	1718.10	1721.8	0614.8	0710S	0710	0712











	FLAT END CYLINDER			ROUND END CYLINDER			MOD. BEVELED CYL.		BEVELED CYLINDER		
											
Order / ISO #	1-837.012	1-837.014	1-837.016	1-881-009	1-881.012	1-881.014	1-878.012	1-878.014	1-885.010	1-885.012	1-886.012
Grits	C	C	C	C F	SC	SC	C F	C	C M	C F	C M
Head / Lgh / Tip (mm)	1.2/7/1.2	1.4/8/1.4	1.6/8/1.6	0.9/7/0.9	1.2/7/1.2	1.4/8/1.4	1.2/8/1.2	1.4/8/1.4	1.0/8/1.0	1.2/8/1.2	1.2/10/1.2
Neo. #	0712.7		0716.8	1210.7	1212.7	1214.8	1800.8	1804.8		1812.8	1812.10

	BEVELED CYLINDER			ROUND				INVERTED CONE				
												
Order / ISO #	1-886.016	1-801.010	1-801.012	1-801.014	1-801.016	1-801.018	1-801.021	1-801.023	1-805.012	1-805.016	1-807.016	1-807.018
Grits	C	M	M	C	C M F	C M	C	C M F	M	C	C	C
Head / Lgh / Tip (mm)	1.6/10/1.6	1.0/0.8	1.2/1.0	1.4/1.2	1.6/1.4	1.8/1.6	2.1/1.9	2.3/2.0	1.2/0.9	1.6/1.4	1.6/4.0	1.8/5.0
Neo. #		0110	0112		0116	0118		0123	0312	0316	0316.4	0318.5

	DOUBLE INVERTED CONE			PEAR				END CUTTER			INTER-PROXIMAL
											
Order / ISO #	1-813.014	1-S830.008	1-830.008	1-830.010	1-S830.012	1-830.012	1-830L.012	1-10839.012	1-10839.014	1-10839.016	1-392.016
Grits	C	M	M	M	M	M	C	M	M	M	M F
Head / Lgh / Tip (mm)	1.4/2.2	0.8/1.6	0.8/1.6	1.0/2.7	1.2/2.7	1.2/2.7	1.2/4	1.2	1.4	1.6	1.6/5
Neo. #	0416	0508S	0508		0512S	0512		5012	5014	5016	1416

	NEEDLES							WHEEL		EGG
										
Order / ISO #	1-S858.008	1-858.008	1-859L.012	1-858.010	1-858.012	1-859.014	1-859.016	1-909.035	1-909.042	1-379.023
Grits	F	M F XF	C F	M	C M F XF	C M F XF	C F XF	C	C	SC M F XF
Head / Lgh / Tip (mm)	0.8/3	0.8/3	1.2/11	1.0/8	1.2/8	1.4/10	1.6/10	3.5/1.3	4.2/1.5	2.3/4.2
Neo. #	1300FS	1310.3	1312.11		1312.8	1314.10		2035	2042	1900

	EGG		FOOTBALL				DEPTH CUTTER			FLAME		
												
Order / ISO #	1-379.018	1-368.016	1-368.020	1-368.023	1-828.026	1-834.016	1-834.021	1-889.009	1-862.010	1-S862.012	1-861.012	1-862.012
Grits	C F	C F XF	C F	SC M F XF	M	M	M	M F	C M F	C M	C XF	C M F XF
Head / Lgh / Tip (mm)	1.8/3.4	1.6/3.2	2.0/5.0	2.3/5.2	dep: 0.5	dep: 0.3	dep: 0.5	0.9/3.0	1.0/8.0	1.2/8.0	1.2/7.0	1.2/8.0
Neo. #	1908	1916	1920	1923	6005	7003	7005	1509.3	1510.8	1512.6	1512.7	1512.8

	FLAME					OCCLUSAL REDUCTION		GROSS REDUCTION		ACORN
										
Order / ISO #	1-863.012	1-861.014	1-862.014	1-862.016	1-862L.018	1-811.033	1-811.037	1-6051.17	1-6055.18	1-905.028
Grits	C M F XF	C M	C	C M F XF	C F	SC	SC	C	C	M F
Head / Lgh / Tip (mm)	1.2/10	1.4/7.5	1.4/8	1.6/8	1.8/9.0	3.3/5.0	3.7/7.0	1.7/8/1.7	1.8/8/1.3	2.8/2.9
Neo. #	1512.10	1514.7	1514.8	1516.8	1518.9	2133	2137	2217	2218	

MULTI-USE OR SINGLE USE | NON-STERILIZED

PENGUIN VALUEDIAMOND ECONOMIC LINE OF DIAMOND BURS

This economic line of diamond burs is very high quality and come in packs of 10 burs.

ValuDiamond burs are designed for both single and multiple use; their one-piece stainless steel shanks undergo a special heat treatment to ensure rigidity, strength and precision.

There are over 200 shapes and grits available, which all sell for the same low price.

10/PACK

REG \$29.99

PROMO BUY 3 PACKS

NET \$23.99/PK

ROUND END TAPER

Order / ISO #	V-850.010	V-856.012	V-850.012	V-S855.014	V-856.014	V-850.014	V-S855.016	V-856.016	V-850.016
Grits	C M	C M	C M	SC C	SC C M F	SC C M F	SC M	SC C M F XF	SC M F
Head / Lgh / Tip (mm)	1.0/10/0.7	1.2/7/0.7	1.2/10/0.7	1.4/6/0.9	1.4/8/0.9	1.4/10/0.9	1.6/6/1.1	1.6/8/1.1	1.6/10/1.1
Neo. #		1112.7	1112.10	1114.6	1114.8	1114.10	1116.6	1116.8	1116.10

ROUND END TAPER

Order / ISO #	V-856.018	V-856N.018	V-856L.018	V-856.021	V-856.022	V-S856.024	V-856.024	V-847.010	V-848.010	V-847.012
Grits	SC C M F	C	SC C M F	SC	SC	SC	SC F	C F	C M	C
Head / Lgh / Tip (mm)	1.8/8/1.4	1.8/8/1.2	1.8/9/1.4	2.1/8/1.6	2.2/8/1.8	2.4/7/1.8	2.4/7/1.8	1.0/8/0.7	1.0/10/0.7	1.2/8/0.8
Neo. #	1118.7	1118.8	1118.9	1121.8	1122.8	2424CS	2424			

FLAT END TAPER

FLAT END TAPER

Order / ISO #	V-847L.014	V-847.014	V-S847.014	V-847.016	V-848.016	V-847.018	V-848.018	V-S846KR.016	V-847KR.016	V-845KR.018	V-847KR.018
Grits	C	SC C M F	SC	SC C M F	SC M	SC M	SC M	SC M	SC C M F	SC C M	SC
Head / Lgh / Tip (mm)	1.4/9/0.7	1.4/8/1.0	1.4/8/1.0	1.6/8/1.2	1.6/10/1.1	1.8/8/1.4	1.8/10/1.4	1.6/6/1.1	1.6/8/1.1	1.8/4/1.4	1.8/8/1.4
Neo. #	0914.9	0914.8	0914.8	0916.8	0916.10	0918.8	0918.10	0816.6S	0816.8	0818.4	0818.8

MODIFIED FLAT END TAPER

MODIFIED FLAT END TAPER

MODIFIED CHAMFER

Order / ISO #	V-845KR.025	V-878K.012	V-879K.012	V-878K.014	V-879K.014	V-S877K.016	V-878K.016	V-878K.018	V-879K.016	V-879K.018	V-878K.021
Grits	C M	C M F	C M	SC C M F	C M	C M	SC C M F	SC F	SC	SC M	C
Head / Lgh / Tip (mm)	2.5/4/2.0	1.2/8/0.8	1.2/10/0.8	1.4/8/1.0	1.4/10/1.1	1.6/6/1.2	1.6/8/1.2	1.8/8/1.4	1.6/10/1.2	1.8/10/1.4	2.1/8/1.6

FLAT END CYLINDER

MOD. FLAT END CYL.

ROUND END CYLINDER

Order / ISO #	V-835.008	V-835.009	V-S835.010	V-835.010	V-835.012	V-835.014	V-837.010	V-837.012	V-837.014	V-837.016	V-837KR.014	V-881.009	V-881.012	V-881.014
Grits	C	C	C M	C M	C	C	C	C	C	C	C F	C F	SC	SC
Head / Lgh / Tip (mm)	0.8/4/0.8	0.9/4/0.9	1.0/4/1.0	1.0/4/1.0	1.2/4/1.2	1.4/4/1.4	1.0/7/71.0	1.2/8/1.2	1.4/8/1.4	1.6/8/1.6	1.4/8/1.4	0.9/7/0.9	1.2/7/1.2	1.4/8/1.4
Neo. #			0710S	0710	0712			0712.7		0716.8	0614.8	1210.7	1212.7	1214.8

MOD. BEVELED CYL.

BEVELED CYLINDER

ROUND

ORDER/ISO #	V-878.012	V-878.014	V-885.010	V-885.012	V-886.012	V-886.016	V-801.009	V-801.010	V-801.012	V-801.014	V-801.016	V-801.018
GRITS	C F	C F	M	C F	C M	C	M	M	C M	C M	C M F	C M
Head/Lgh/Tip(mm)	1.2/8/1.2	1.4/8/1.4	1.0/8/1.0	1.2/8/1.2	1.2/10/1.2	1.6/10/1.6	0.9/0.8	1.0/0.8	1.2/1.0	1.4/1.0	1.6/1.4	1.8/1.6
Neo. #	1800.8	1804.8		1812.8	1812.10			0110	0112		0116	0118

ROUND

INVERTED CONE

DOUBEL INVERTED CONE

PEAR

ORDER/ISO #	V-801.021	V-801.023	V-805.012	V-805.016	V-807.016	V-807.018	V-813.014	V-S830.008	V-830.008	V-830L.008	V-830.010	V-S830.012
GRITS	C	C M F	M	C	C	C	C	M	M	M	M	M
Head/Lgh/Tip(mm)	2.3/2.0	2.3/2.0	1.2/9	1.6/1.4	1.6/4.0	1.8/5.0	1.4/2.2	0.8/1.6	0.8/1.6	0.8/2.7	1.0/2.7	1.2/2.7
Neo. #	0123	0123	0312	0316	0316.4	0318.5	0416	0508S	0508			0512S

PEAR

END CUTTER

INTER-PROXIMAL

NEEDLES

ORDER/ISO #	V-830.012	V-830L.012	V-10839.012	V-10839.014	V-10839.016	V-392.016	V-S858.008	V-858.008	V-858.010	V-859.010
GRITS	M C	C	M	M	M	M F	F	M F XF	M F XF	M F XF
Head/Lgh/Tip(mm)	1.2/2.7	1.2/4	1.2	1.4	1.6	1.6/5	0.8/3	0.8/3	1.0/8.0	1.0/10
Neo. #			5012	5014	5016	1416	1300FS	1310.3		

NEEDLES

0512S

0512

WHEEL

EGG

5016

ORDER/ISO #	V-858.012	V-859.012	V-858.014	V-859.014	V-859.016	V-909.035	V-909.042	V-379.023	V-379.018
GRITS	C M F XF	C F XF	C M F XF	C M F XF	C M F XF	C	C	SC M F XF	C F
Head/Lgh/Tip(mm)	1.2/8.0	1.2/11	1.4/8	1.4/10	1.6/10	3.5/1.3	4.2/1.5	2.3/4.2	1.8/3.4
Neo. #		1312.11	1314.8	1314.10		2035	2042	1900	1908

	FOOTBALL				DEPTH CUTTER			FLAME			
ORDER/ISO #	V-368.016	V-368.018	V-368.020	V-368.023	V-828.026	V-834.016	V-834.021	V-889.009	V-862.010	V-S862.012	V-861.012
GRITS	C F XF	C	C F SC	C M F XF	M	M	M	M F	C M F XF	C M	C F XF
Head/Lgh/Tip(mm)	1.6/3.2	1.8/4.0	2.0/5.0	2.3/5.2	dep: 0.5	dep: 0.3	dep: 0.5	0.9/3.0	1.0/8.0	1.2/8.0	1.2/7.0
Neo. #	1916		1920	1923	6005	7003	7005	1509.3	1510.8	1512.6	1512.7

	FLAME						OCCLUSAL REDUCTION		GROSS REDUCTION		ACORN
ORDER/ISO #	V-862.012	V-863.012	V-861.014	V-862.014	V-862.016	V-862L.018	V-811.033	V-811.037	V-6051.17	V-6055.18	V-905/028
GRITS	C M F XF	C M F XF	C M	C M F XF	C M F XF	C F	SC	SC	C	C	M F
Head/Lgh/Tip(mm)	1,2/8.0	1.2/10	1.4/7.5	1.4/8	1.6/8	1.8/9.0	3.3/5.0	3.7/7.0	1.7/8/1.7	1.8/8/1.3	2.8/2.9
Neo. #	1512.8	1512.10	1514.7	1514.8	1516.8	1518.9	2133	2137	2217	2218	

NON STERILE

PENGUIN CARBIDE

PROMO BUY 3 +

FG CARBIDES

10 PER PKG. \$29.99 - NET 24.9

PH-1-FG2	PH-1-FG4	PH1-FG	PH1-FG8	PH1-FG245	PH1-FG330
PH-1-FG331	PH-1-FG557	PH-1-FG1157	PH1-FG1557	PH-1-FG1558	

RA CARBIDES

10 PER PKG. \$29.99 - NET 24.95

FGOS CARBIDES

10 PER PAK 59.99

PH-1-RA2	PH-1-RA4	PH-1-RA6	PH1-RA8	PH-1-FGOS557	PH-1-FGOS702

ETHICON SUTURES



42497	VICRYL FS-2, 4-0 UNDYED 27", 36/BX	J422H	\$269.99
46843	VICRYL 4-0, FS-2, 18", 36/BX	J392H	\$269.99
12471	VICRYL 4-0, P-3, 3/8 CIRCLE, 18", 12/BX	VR494	\$179.99
12470	VICRYL 3-0, FS-2, 3/8 CIRCLE, 27", 36/BX	VR2298	\$319.99
14061	VICRYL 4-0, PS-2, 18", 12/BX	J496G	\$159.99
14044	VICRYL 3-0, PS-2, 18", 12/BX	VR497	\$159.99
12472	VICRYL RAPIDE 3-0, PS-2, 18", 12/BX	632G	\$179.99
46841	SILK-XI NEEDLE 18" SILK, 12/BX	612G	\$69.99
40676	SUTURE CUTICULAR X-1 SZ 3 18" GUT PLAIN, 12/BX	635H	\$104.99
28834	ETHICON GUT FS-2, 4-0 CHROMIC 27", 36/BX	683G	\$299.99
42230	ETHICON GUT FS-2, 4-0 CHROMIC 27", 36/BX	1654G	\$69.99
28835	CHROMIC GUT 4-0, P-3 NEEDLE, 18", 12/BX	687G	\$179.99
28840	CHROMIC GUT 5-0, P-3 NEEDLE, 18", 36/BX	1654G	\$194.99
12400	PLAIN GUT 4-0, FS-2 NEEDLE, 21", 36/BX	687G	\$299.99
12772	PROLENE BLUE 4-0, P-3 NEEDLE, 18", 12/BX	H821H	\$179.99
24089	PROLENE BLUE 5-0, P-3 NEEDLE, 18", 12/BX	8698G	\$179.99
24090	PROLENE BLUE 6-0, P-3 NEEDLE, 18", 12/BX	8695G	\$179.99
24092	PROLENE BLUE 6-0, P-1 NEEDLE, 12", 12/BX	8697G	\$179.99
12495	ETHILON FS-2, 5-0, 18", 12/BX	6616G	\$69.99
10871	SILK FS-1, 3-0, 18", 12/BX	684G	\$69.99
43787	VICRYL P-3, 5-0, 18", 12/BX	J493G	\$159.99
50031	SILK C-3, 4-0, 27", 12/BX	736G	\$94.99

SURGICEL ABSORBABLE HEMOSTAT

39526
5.0 CM X 1.2 CM
(JJ-1955)



\$199.99
BOX OF 12

ABSORBABLE HEMOSTAT ABSORBABLE HEMOSTAT GELATIN SPONGE

1 CM CUBE, 24/BOX
PH-MS0005

~~\$69.99~~
PROMO \$55.99

GELATIN SPONGE ABSORBABLE GELATIN SPONGE


32262
BOX OF 12 ENVELOPES - STERILE
12-7, 2 CM X 6 CM
(12 CM SQ.) X 7 MM (THICK)
(#1972)



\$139.99
BOX OF 12

ETHICON

4+1



ALM SCALPEL BLADES



DISPOSABLE SCALPEL S/S
AVAILABLE IN #11, 12, 15 (10/BX) **\$12.99**

SURGICAL BLADES S/S
AVAILABLE IN #11, 12, 15 (100/BX) **\$29.99**

SURGICAL BLADES S/S
15C (100/BX) **\$45.99**

OTHER SIZE BLADES AVAILABLE, PLEASE CALL. FINAL SALE. NO RETURNS.

JASPER™ ISOLATION DISPOSABLE YELLOW GOWNS

PACK OF 10



REORDER CODE
PH-IDG-Y10

NET \$7.79 /PK OF 10

REG **\$12.99 /PK OF 10**

LOUISE™ LEVEL 2 DISPOSABLE GOWN

PACK OF 100



PROMO

NET **\$156 /CASE**

PROMO

REG **\$195 /CASE OF 100**

REORDER CODE
PH-L2100

MORAINE™ LEVEL 3 WASHABLE GOWNS

8+8



NET **\$7.47 /PACK**

8+8

REG **\$14.95 /PACK**

REORDER CODE
L31502

EVEREST™ LEVEL 3 WASHABLE GOWN

8+8



NET **\$17.24 /PACK**

8+8

REG **\$34.49 /PACK**

REORDER CODE
L3150

JASPER™ SHOE COVERS

10% DISCOUNT



NET **\$17.99 /BOX**

10% DISCOUNT

REG **\$19.99 /CASE OF 100**

REORDER CODE
PH-SC100

JASPER™ BOUFFANT CAP

10% DISCOUNT



NET **\$17.55 /BOX**

10% DISCOUNT

REG **\$19.5 /CASE OF 100**

REORDER CODE
BC100

JASPER™ REUSABLE FACESHIELD

10% DISCOUNT



NET **\$17.55 /BOX**

10% DISCOUNT

REG **\$19.5 /CASE OF 100**

REORDER CODE
FS-PK10

SAFETY GOGGLES SINGLE UNIT

10% DISCOUNT



NET **\$17.55 /BOX**

10% DISCOUNT

REG **\$19.5 /UNIT**

REORDER CODE
SG0077

HB FACIAL TISSUE 30X100/CS

10% DISCOUNT



NET **\$44.99 /BOX**

10% DISCOUNT

REG **\$49.99 /1 BOX**

REORDER CODE
FT10030

CASCADES PRO SELECT™ WHITE MULTIFOLD PAPER TOWEL, 16X250 SHEETS

10% DISCOUNT



NET **\$59.95 /BOX**

10% DISCOUNT

REG **\$64.99 /1 BOX**

REORDER CODE
MF4000K

CASCADES PRO SELECT™ KRAFT MULTIFOLD HAND PAPER TOWEL, 16X250 SHEETS

10% DISCOUNT



NET **\$49.95 /BOX**

10% DISCOUNT

REG **\$54.99 /1 BOX**

REORDER CODE
PAP108

HiTeck.

Your Partner in Dentistry


ORAL SURGERY INSTRUMENTS



SCAN TO SEE ALL THE
INSTRUMENT OPTIONS

BUY 5 GET 1 FREE

Free Instrument must be of same
or lesser value



HYGIENE & EXAM INSTRUMENTS

**BUY 5 GET
1 FREE**

Free Instrument must be of same or lesser value

PRODUCT NAME	ITEM CODE	PRICE
204S Sickle Scaler	1184	34.99
H6/H7 Hygienist Scaler	1175	34.99
Super Taha Scaler	STAHA	34.99
Taha 1 Anterior Sickle Scaler	TAHA1	34.99
Taha 2 Posterior Sickle Scaler	TAHA2	34.99
Taha 3 Posterior Sickle Scaler	TAHA3	34.99
Taha 4 Posterior Sickle Scaler	TAHA4	34.99
Super Taha Scaler	STAHA	34.99
11/12 Gracey Curette	1234	34.99
13/14 Gracey Curette	1235	34.99
15/16 Gracey Curette	1236	34.99
1/2 Gracey Curette	1229	34.99
1/2 Mini Rigid Gracey Curette	G12MR	34.99
1/2 Arhum Gracey Curette	AG12	34.99
11/14 Arhum Gracey Curette	AG1114	34.99
13/14 Rigid Gracey Curette	G1314R	34.99
15/16 Rigid Gracey Curette	G1516R	34.99
7/8 Younger Good Universal Curette	1212	34.99
4R/4L Columbia Universal Curette	1216	34.99
13/14 Columbia Universal Curette	1217	34.99
H5/J33 Hygienist Jacquette Sickle Scaler	1176	34.99
5/6 Langer Universal Curette	1206	34.99
1/2 Langer Universal Curette	1204	34.99
13/14 Indiana University Universal Curette	1208	34.99
ODU 11/12 Elongated Explorer	1122AF	29.99
ODU 11/12 Explorer	1122	29.99
Signature Mirror Handle, Cone Socket	1168	21.99
College Utility Dressing Plier	1157	19.99
Meriam Utility Dressing Plier	1158	19.99
23/UNC15 Expro	1151	29.99
23/QW Williams Expro	1146	29.99
23/CP 12 Siyah Lachak Expro	1150S	39.99
UNC15/11.5B Probe	1137	29.99
18 Ball Burnisher, 1.6MM 2.4MM	1918	32.99
1 Goldstein Plastic Filling Instrument, Siyah Series	2467	42.99
2 Goldstein Plastic Filling Instrument, Siyah Series	2468	42.99
3 Goldstein Extra Flex Plastic Filling Instrument, Siyah Series	2469	42.99
4 Goldstein Extra Flex Plastic Filling Instrument, Siyah Series	2470	42.99
27/29 Ball Burnisher, 1.6MM 3.1MM, Siyah Series	300-104	42.99
Small/Medium Ball Burnisher, Siyah Series	300-132	44.99
T-Flex Composite Instrument	TFLEX	98.99
E2 Spoon, 1.5MM Excavator	1808	32.99
#18 Spoon, 1.5MM Excavator	1813	32.99
26/27S Ball Burnisher, 1.3MM 2.1MM	1919	32.99
28/29 Ball Burnisher	1921	32.99
L7 Levy Carver	2397	32.99
1/2 Hollenback Carver, Siyah Series	1904S	42.99
Interproximal Carver, Long, Siyah Series	2369	42.99

RAYSCAN $\alpha+$ 8x8 12x10
16x10

3D CBCT



FINANCING FOR
36 MONTHS

All-In One

A Digital Panoramic
X-Ray, cone Beam CT,
Cephalometric X-Ray,
and Impression Scan.

\$2490/month



RAYSCAN $\alpha+$ 2D



FINANCING FOR
36 MONTHS



All - In One

A Digital Panoramic
X-Ray

\$999/month

Comfortek® Dental Operating System

Comfortek X1 Series dental operating system is tailored for professionals, providing more efficiency during treatments and improving the working environment.



Reliable & Durable

Main frame is composed of cast aluminum and steel parts.



Main Switch

Reduces the chances for accidents during off-time with a main switch for all (air, water and electricity).



Flexible Control

The wide range swivel cuspidor bowl and assistant element provides a more spacious treatment area.



Hygienic

Cleaning and disinfection with easy removable suction filter kit & 360° rotatable stainless water outlet.



Comfortable and Considerate

With double-articulated headrest you can easily adjust for comfort.



Warranty

1 Year Full Warranty
5 Years on motors

Color Option:



Espresso Brown



Carbon Black



Burgundy



Cobalt Blue



Sky Blue



Medium grey



Off-white



Olive Green



Mustard Yellow



Comfortek X1® EXPERT

Dental Operating System

INCLUDES

- Stellar patient chair
- Luxury upholstery
- Delivery unit
- Cuspidor
- Floor standing water unit
- L&R long arm rest
- Assistant and doctor control panel with memory position
- M5 LED light
- Built in LED piezo scaler
- Premium doctor stool and assistant stool
- 3 year warranty on parts
- 7 year warranty on motor



FINANCING FOR 36 MONTHS

Comfortek X1 EXPERT
(20401)

\$469/MONTH

Comfortek X1®

Comprehensive & Advanced Technology

INCLUDES

- Stellar patient chair
- Luxury sewn upholstery
- Deluxe Delivery unit
- Floor standing water unit
- L&R long arm rest
- Cuspidor
- Assistant and doctor control panel with memory position
- M8 LED operating light, white and yellow light sensor
- 2 Fiber-optic tubing
- Added electric motor
- Built in LED piezo scaler
- Premium doctor stool and assistant stool
- 3 year warranty
- 7 year warranty on motor



FINANCING FOR 36 MONTHS

Comfortek X1 MASTER
(20504)

\$630/MONTH

Optional Configurations:



Proportional Foot Control



Air Lock Switch with 4 Way Joystick



Second HVE Suction



Monitor Arm



Quick Disconnect (Air or Water) x1



Quick Disconnects (Air and Water) x2

and more...



30-Day
Satisfaction
Guarantee



I-SENSOR

Intraoral Sensor

Available in 3 sizes



(20231) 31*21*4.4mm
Reg. \$3950
Sale: \$3,250



(19036) 34*25*4.4mm
Reg. \$3950
Sale: \$3,510



(20230) 37*27*4.4mm
Reg. \$4250
Sale: \$4,095



(20232)

BUNDLEPACKAGE BUY ANY 2 SENSORS

Get \$500 Gift Cards or \$700 Penguin Credit

NET \$6888

2+1

REG \$8200

FLIGHT A12 OPERATORY CHAIR PACKAGE

(INCLUDES FREE DENTAL STOOL)

\$599/month

0%
FINANCING FOR
48 MONTHS



Features:

- A12 Hydraulic Patient Chair
- State of the Art Hands free Sensor Controlled Flight Torch LED Light
- Synchronized Chair and Back Movement for superior Ergonomics
- Choose from 9 Plush Sewn or Seamless Upholstery Colors and 2 Ultraleather Colors
- Self Contained Water Bottle System
- Thin Backrest and Streamlined Seat rest
- Designed and Engineered in Canada
- DCI Compatible Parts

AUTO CLAVE B

\$330/month

The **Clave B** Autoclave features a pre/post vacuum that provides full penetration of steam into wrapped, porous, and hollow loads. The pre-vacuum phase enables the complete removal of air from the chamber while the post-vacuum air removal ensures that the load is completely dry.

The combination of the pre/post vacuum allows closed door drying and sterilization for the most difficult loads.

The Clave B is an intuitive and easy to use Class B Sterilizer. All operations are controlled with only 4 buttons, with 3 buttons that control the program selection and 1 to start the sterilization cycle.

Health Canada and ETL/CSA Approved
Excellent Sterilizer for use in Dental Offices, Veterinarian Offices, Medical Offices, Podiatry, Chiropractist, Tattoo and Beauty Salons

0%
FINANCING FOR
36 MONTHS



PROFESSIONAL EXPERTS IN DENTAL CONSTRUCTION

BRINGING YOUR VISION TO LIFE

[REQUEST A QUOTE](#)



DENTAL CONSTRUCTION – FUNCTIONAL AND BEAUTIFUL DENTAL OFFICE CONSTRUCTION

Penguin Dental Construction is your trusted partner in dental office construction, renovation, and remodeling in Ontario. With our expertise in the field and a commitment to excellence, we are dedicated to creating functional and aesthetically pleasing dental environments. Whether you're starting from scratch or looking to upgrade your existing clinic, our specialized design and construction services cater specifically to the unique needs of the dental industry.

[REQUEST A QUOTE](#)

OUR SERVICES

- NEW BUILDS
- RENOVATION
- IPAC CONSULTING ONTARIO
- DENTAL X-RAYS ROOM DESIGN
- DENTAL OFFICE DESIGN

OUR PROCESS

- BUILDING SOLUTIONS
- PROJECT MANAGEMENT
- VALUE-ADDED APPROACH

CONSTRUCTION PROMO

FREE CONSULTATION, FREE OFFICE DESIGN,
AND DISCOUNT ON EQUIPMENT!



PENGUIN
HEALTH