



Q4 2026 PROMOTION CATALOGUE

OFFERS VALID TO DEC 31, 2026

PENGUIN HEALTH PRESENTS

ATLAS™

Elevating your sterilization to new heights.

Sterilization, preventative, restorative, disposables and equipment — this quarter's best offers from one Canadian supplier.

SUMMIT BUNDLE

In-clinic exclusive savings

20-MINUTE RESULTS

ATLAS Rapid Incubator

FREE IPAC COURSES

For Summit Bundle clients

PROUDLY CANADIAN

Health Canada registered



A PROUD CANADIAN COMPANY penguinhealth.ca • 905-901-1579 • orders@penguinhealthgroup.com

Contents

03

Who we are

Penguin Health — proudly Canadian

04

★ SPECIAL PROMO: FREE IPAC Course & Sterilization Bundle

The ATLAS Summit Bundle — save up to \$5,256.74

05–15

ATLAS Sterilization

Rapid incubator, indicators, ATLAS Wash, Tracker, packaging, cleaning, equipment, IPAC

16–17

Preventative — LUNA Whitening

NEW LUNA Practice Starter Bundle \$399

18–19

Restorative

Precision Edge · O-Star curing light \$959

20–28

Harvard — Made in Germany

MultiChrome, UltraFill, BioCal-CAP, InterLock, ZirconCore

29–30

Penguin Sterilization

Pouches, wraps, barrier film, surface & hand hygiene

31–32

Penguin Daily Disposables

Gloves, tips, ejectors, bibs, prophyl cups, applicators

33–34

Preventative — Woodpecker

PT-A air polisher, scaler tips and powders

35

Penguin Handpieces

Low-speed, prophyl and air motors

36–37

HiTeck Instruments

Oral surgery and hygiene & exam instruments

38–40

Capital Equipment

Chairs, sterilizers, Smile 2 scanner \$14,999

41–43

Office Partners, Notes & Order Form

Turnkey construction promo and quick order form

NET \$

4+1

REG \$

PACK

How to read our prices: NET = your price per unit with the full offer. REG = regular price. PACK = what one unit contains. Prices in CAD, before taxes & freight.



PENGUIN HEALTH • PROUDLY CANADIAN

Canadian-made sterilization monitoring.

One Canadian supplier for every link in your sterilization chain — and one team to call.

Canadian

Headquartered in
Mississauga, Ontario

One supplier

Indicators,
instruments, records
& training

Health Canada

Registered products

Lifetime

Warranty on all our
products



Full documentation

IFUs, SDS and specifications
for every product



Lot traceability

Lot & expiry printed on every
indicator



Training & support

From the people who specify
the product



ATLAS Rapid



Biological indicators



Chemical indicators



Cleaning verification



Digital records

★ SPECIAL PROMO

IN-CLINIC EXCLUSIVE

The ATLAS Summit Bundle

FREE

\$2,000 VALUE

INCLUDED WITH THE SUMMIT BUNDLE

FREE IPAC Sterilization Assessment

2-hour on-site review by a CIC-certified consultant, with a written report. Or \$2,000 off the 3.5-hour Assessment & Training — just \$1,500.

WHAT'S INCLUDED · CODE	QUANTITY & PACK SIZE	IN THE BUNDLE
 ATLAS Rapid Incubator PH01980	1 incubator	FREE*
 ATLAS Ultra Rapid BI — 20 min AX224	4 boxes · 50 per box · 200 spore tests	✓ INCLUDED
 ATLAS Shield Type 5 Integrator AT-AXCP90	20 packs · 250 per pack · 5,000 strips	✓ INCLUDED
 ATLAS Wash-U-Test PH-AT-50	4 packs · 50 per pack · 200 tests	✓ INCLUDED
 ATLAS Shield PCD — stainless steel PH-2026	1 device	FREE
 Wash-U-Test Stainless Holder PH-AT-01H	1 holder	FREE
 IPAC Readiness Assessment 2-hour on-site visit · CIC-certified	1 visit for your team	FREE

REGULAR PRICE

~~\$7,248.70~~

Products + \$2,000 IPAC assessment

SUMMIT BUNDLE PRICE

\$1,991.96

Quote SUMMIT BUNDLE

YOU SAVE

\$5,256.74

72% off the regular price

\$3,256.74 off the products + FREE \$2,000 IPAC assessment

TOTAL SAVINGS \$5,256.74

*Incubator free with trade-in of a qualifying biological incubator, or under the buy-back programme. Limited time, in-clinic only. Validated for genuine ATLAS indicators only.

Book your visit: 905-901-1579

FIVE LINES · ONE CHAIN OF ASSURANCE

Every link in reprocessing, monitored and documented by **one system.**

01



Clean

Verify the ultrasonic and washer are removing soil.

ATLAS Wash

02



Pack

An integrator inside every package.

Shield Type 5

03



Sterilize

Daily air-removal test for Class B autoclaves.

Shield Bowie-Dick

04



Monitor

Prove lethality with a resistant spore population — result in 20 minutes.

ATLAS BIs + Rapid

05



Document

Digitally Track each load to the patients it served.

ATLAS Tracker

One supplier. One conversation when it matters.

Indicators, cleaning verification, packaging and digital records from one Canadian supplier — with **one team to call.**

Standards & documentation

Type 5 integrators comply with **ISO 11140-1** and carry cycle references, lot and expiry. IFUs ship with every instrument; SDS and registration documentation on request.



Follow every step to meet Canadian health standards (CSA Z314).



A proud Canadian company · Health Canada registered products
info@penguinhealthgroup.com · +1 905-901-1579 · penguinhealth.ca

MINI AUTO ULTRA RAPID BIOLOGICAL INCUBATOR

ATLAS Rapid Incubator

PH01980

AXRM100

20
MINUTES


THE COMPLETED PRINTOUT

ATLAS RAPID
BI INCUBATION REPORT
PH01980 / AXRM100

Date	2026-09-30
Well	02
Program	STEAM 60°C
BI code	AX224
BI lot	A1000
Sterilizer	STERI-01
Cycle no.	0142
Start	09:02
End	09:22

RESULT: NEGATIVE

PASS · LOAD RELEASED

Signature

★ **BUILT-IN THERMAL PRINTER**
Every result printed & time-stamped

Sample printout, for illustration

A compact benchtop auto-readout incubator with **four independently timed wells**, a 3.5" colour touch screen and a built-in thermal printer. Every result is **time-stamped, printed and stored** against the sterilizer and operator that produced it.

4

Independent wells

3.5"

Colour touch screen

37/60°CSelectable
temperature**10-15
min**To reach set
temperature**Lifetime**Warranty on all
products

HOW IT READS

1 · Activate

Crush the BI ampoule in the built-in crusher after the cycle.

2 · Insert

Place it in any free well; each well starts its own countdown.

3 · Read

The indicator is read optically — no colour judgement by staff.

4 · Print

Result printed and stored; export over USB for your records.





HERO PRODUCT · SELF-CONTAINED

ATLAS AX224 Ultra Rapid Steam BI

AX224

G. stearothermophilus

Box of 50

The rapid-readout indicator the ATLAS Rapid Incubator is built to read — **20 minutes at 60°C**. A calibrated spore population and its growth medium sit in one sealed unit with a crushable ampoule: nothing to prepare, measure or store, and no opportunity to contaminate the culture in handling.

- **20-minute result** — release loads in the same shift, not tomorrow
- **Self-contained** — spores and medium in one vial, nothing to prepare
- **Lot & expiry** on every vial; results printed and filed automatically

HOW IT WORKS

01

Place

Put the AX224 in a test pack or the hardest-to-reach spot in the load. Label it with the load and date.

02

Sterilize

Run the steam cycle as normal. Steam must reach and kill the spores inside the vial.

03

Activate

After the cycle, let it cool, then crush the ampoule so the growth medium soaks the spore strip. Activate one unprocessed vial from the same lot as a control.

04

Read in 20 min

Place both in the ATLAS Rapid Incubator at 60°C. The result reads automatically in 20 minutes and is printed and stored.

WHAT IT IS

A **biological indicator (BI)** is the only test that proves sterilization actually killed living organisms. Each AX224 vial holds a calibrated population of *G. stearothermophilus* spores — the organism most resistant to steam — with its growth medium sealed in a crushable glass ampoule.

Pass: no growth — the spores were killed. **Fail:** growth detected — quarantine the load and do not use the sterilizer until it is retested.

WHEN TO USE IT

- At least daily, in each sterilizer that is in use
- In every load containing implants — release only after a pass
- After installation, relocation or major repair of the sterilizer
- After any failed BI, until the sterilizer passes again

Confirm frequency with your IPAC team and CSA Z314.

ALSO
AVAILABLE



STEAM · AT-1092

ATLAS Shield™ BI — 10 hr

For sites incubating overnight.



STEAM · AT-2492

ATLAS Shield™ BI — 24 hr

Standard 24-hour steam BI.



Three routes. *One summit.*

Pick the readout speed that fits how your practice already works.

ATLAS RAPID • MODEL PH01980

20 minutes changes everything

Incubated, read and time-stamped for you. Spore results in 20 minutes instead of overnight — implant and critical loads released in the same appointment.

4 INDEPENDENT WELLS

THERMAL PRINTER

USB EXPORT

LIST PRICE

~~\$2,049.00~~

WITH BUY-BACK

\$1,499.00

SUMMIT BUNDLE

FREE*



09:00

Cycle ends

The load comes out with an ATLAS indicator in it.

09:02

Activate

Crush the ampoule and place it in the ATLAS Rapid.

09:22

Result

Read and time-stamped automatically.

09:25

Released

Same appointment, no overnight quarantine.

**4+1
BUY 4 GET 1
FREE**



ATLAS AX224 — Ultra Rapid Steam BI

The rapid-read, self-contained indicator for the ATLAS Rapid. Spore population and growth medium in one crushable ampoule — nothing to prepare or store. *G. steartothermophilus* • box of 50.

NET **\$311.99** / box of 50

REG \$389.99

SAVE \$78.00

Bundle of 5 boxes (250 spore tests): **\$1,559.96**

REORDER AX224

**2+1
BUY 2 GET 1
FREE**

AT
2492



ATLAS Shield™ BI — 24 Hour

Conventional overnight cycle, same construction and lot documentation. Box of 100.

NET **\$199.99** / box of 100

REG \$299.99

SAVE \$100.00

REORDER AT-2492



ATLAS Shield™ BI — 10 Hour

Run at day's end, result next morning. Universal. Box of 50.

PRICE **\$199.99** / box of 50

REORDER AT-1092

*FREE with the Summit Bundle and qualifying trade-in. Buy-back price \$1,499. Validated for genuine ATLAS indicators only. Prices in CAD.

Chemical indicators & air-removal testing

Per-package and per-load evidence between biological tests.

HERO PRODUCT · TYPE 5

AT-AXCP90

 NET
\$26.99

12+8

 REG
\$44.99

PACK OF 250 STRIPS



ATLAS Shield™ Type 5 Integrating Indicator — 5S

Reacts to time, temperature and saturated steam together — place one inside every package. Complies with ISO 11140-1.

59 × 15 mm · pack of 250

134°C 3.5 min · 121°C 15 min

Lot & expiry on every strip



PACK OF 250 · PH-3421A

Penguin Type 5 Moving Integrator

A chemical pellet melts and travels along a wick during the cycle. **Pass when the dark bar reaches the PASS window** — read by position, not shade, so every staff member reads it the same way.

PACK OF 250 · PH-3421A

 NET
\$59.99

4+1

 REG
\$74.99


BOX OF 100 · AX-BD100

ATLAS Shield™ Bowie-Dick Test Cards

FOR CLASS B (PRE-VACUUM) AUTOCLAVES

- **When:** daily, first cycle, empty chamber — before any load.
- **Run:** a Bowie-Dick cycle at 134°C for 3.5 min.
- **Pass:** the pattern turns from pink to an even, uniform black.
- **Fail:** a lighter patch, usually in the centre — air was left in. Stop, retest and call service.
- Record the result on the back · ISO 11140-1 & 11140-4 Type 2.



Nothing sterilizes through debris.

The first link in the chain — and the one **almost nobody tests**

HERO PRODUCT · ULTRASONIC

PH-AT-50

ATLAS WASH-U-TEST®
Cleaning Indicator
Lot: 123456-01
Expiry Date: 11-2023

OK

Medical Care ISO/TS 16949:16

3+1

NET

\$112.49

REG \$149.99

PACK OF 50 PH-AT-50

ATLAS Wash-U-Test™ Ultrasonic Indicator

A coloured test coating is removed only when cavitation is genuinely occurring across the tank — catching a failing transducer, low fill, exhausted solution or a poorly positioned basket. Pack of 50.

Confirms cavitation is happening

Catches transducer failure

Quick weekly check, dated record

ATLAS AMERISCO
Cleaning Indicator
Lot number
Expiry Datm/yyyy

OK

MADE IN USA

PACK OF 250

ATLAS Washer Test Indicator

Complete soil removal confirms the washer's action, chemistry and temperature.

PACK OF 250 · PH-ATWR

NET **\$112.49**

3+1

REG **\$149.99**

STAINLESS STEEL · PH-AT-01H

Wash-U-Test Indicator Holder

Holds the indicator in the same tank position every time, so results compare week to week.

EACH · PH-AT-01H

Q4 2026 PRICE **\$79.99**

PROCESS CHALLENGE DEVICE · PH-2026

ATLAS Shield™ PCD

Reusable challenge device that holds an indicator inside a hard-to-reach path to test steam penetration.

EACH · PH-2026

Q4 2026 PRICE **\$59.99**



Make your sterilization process traceable.

ATLAS Tracker records every critical step and integrates directly into the patient file. A touchscreen workstation, label printers and scanners replace the paper binder — and **every report type is included.**



All-in-one computer

Touchscreen hub for every cycle



Label printers

Load label at the point of packaging



Scanners

Capture at point of use, linked to the patient

EVERY REPORT TYPE INCLUDED — NO ADD-ONS

- | | | |
|----------------------------------|--------------------------------|-------------------------------|
| ✓ Sterilization cycle reports | ✓ Biological indicator reports | ✓ Chemical indicator reports |
| ✓ Patient traceability reports | ✓ Washer & ultrasonic logs | ✓ Waterline testing reports |
| ✓ Failure & corrective action | ✓ Staff training & competency | ✓ Equipment maintenance logs |
| ✓ Audit-ready compliance summary | ✓ Custom & scheduled reports | ✓ One-click regulatory export |

Our team provides full set up and tech support.



A proud Canadian company · Health Canada registered products
info@penguinhealthgroup.com · +1 905-901-1579 · penguinhealth.ca

Packaging & barriers

Pouches, wrap, film and tape — everything that holds a sterile load together.

CLASS 4 · 200 / BOX

Sterilization Pouches — Class 4 Self-Sealing

Peel, fold and seal — no heat sealer needed.

Use: single instruments and small sets.

- ✓ Built-in Class 4 process indicator
- ✓ 4 sizes, 2.75×10" to 5.25×11"



7+3

NET

from \$8.39

REG from \$11.99

200 / BOX · 4 SIZES

SMS WRAP · 500 / BOX

HB SMS Wrap

A barrier layer between two strength layers.

Use: surgical, implant and heavy instrument trays.

- ✓ Resists tearing on sharp tray corners
- ✓ 12, 15, 20 & 24" sheets



10% OFF

NET

from \$71.95

REG from \$79.95

500 / BOX · 4 SIZES

CSR WRAP · 500 / BOX

HB CSR Wrap

Value CSR wrap for high tray volumes.

Use: everyday cassettes and routine trays.

- ✓ Lower cost per tray on routine loads
- ✓ 15, 20 & 24" sheets



10% OFF

NET

from \$80.99

REG from \$89.99

500 / BOX · 3 SIZES

SURFACE BARRIER · PH-BF02C

HB Barrier Film Roll

Perforated adhesive film, clear or blue.

Use: light handles, switches and keyboards.

- ✓ Replace between patients — faster than disinfecting
- ✓ Perforated sheets tear off cleanly



3+1

NET

\$24.74

REG \$32.99

1,200 SHEETS · PH-BF02C

EXTERNAL PROCESS INDICATOR

ATLAS Autoclave Tapes

Seal and mark wrapped packs in one step.

Use: closing SMS/CSR-wrapped trays and packs.

- ✓ Changes colour on exposure to steam
- ✓ Processed and unprocessed packs told apart at a glance



Cleaning & disinfection

From the ultrasonic bath to the operatory surface — from one supplier.

1 US GAL · PH-100-0092

Mark3 General Purpose Ultrasonic Solution

Non-ammoniated liquid for the ultrasonic bath.

Use: daily ultrasonic bath changes.

- ✓ Non-ammoniated — lifts debris from hinges
- ✓ Verify the bath with ATLAS Wash-U-Test



BUY 4+

NET
\$39.99

REG \$49.99

1 GALLON (4 L) BOTTLE

64 TABS · M3100-7640

Mark3 Bioenzymatic Ultrasonic Cleaner

Pre-measured enzymatic tablets in sachets.

Use: heavy protein soil like blood.

- ✓ Same dilution every time — no measuring
- ✓ Anti-corrosive formula



5+1

NET
\$58.32

REG \$69.99

BOX OF 64 TABLETS

1-MIN CONTACT · PH-OS1CA-W12

Optim® OS1 Wipes

Ready-to-use disinfectant cleaner wipe.

Use: fast turnover between patients.

- ✓ 1-minute contact time
- ✓ Clean and disinfect in one step



BUY 24+

NET
\$14.99

REG \$19.99

CAN OF 160 WIPES

160 / CAN · PH-SCI33W12

Optim® 33TB Wipes

Surface cleaner & intermediate-level disinfectant.

Use: surfaces needing a TB claim.

- ✓ 1-min tuberculocide, virucide & bactericide
- ✓ DIN 02285347 · 160 wipes



BUY 24+

NET
\$16.49

REG \$19.95

CAN OF 160 WIPES



FLIGHT DENTAL SYSTEMS

Clean it. Sterilize it.

The two machines at the heart of the sterilization room — monitored end to end with ATLAS.

0%
FINANCING
AVAILABLE



Q4 pricing: ask your representative

STEP 1 · CLEAN · UC-6065

Flight Ultrasonic Cleaner

A heated stainless-steel ultrasonic bath that loosens debris from instruments before packaging — reaching hinges, serrations and lumens a brush cannot.

When to use it: first step after use — before instruments are inspected, packaged and sterilized. Run daily with Mark3 solution or tablets.

KEY FEATURES & BENEFITS

- Heated tank with temperature dial — warm solution cleans faster
- Timer dial for repeatable cycles
- Degas mode removes dissolved air for stronger cavitation
- Sweep mode spreads energy evenly across the tank
- Stainless-steel tank with lid and carry handles
- Verify every bath with ATLAS Wash-U-Test



REG
\$11,769.50

\$330/month
0% FINANCING · 36 MONTHS

STEP 2 · STERILIZE · CLASS B

Flight Clave B Autoclave

A Class B pre- and post-vacuum steam sterilizer for wrapped, porous and hollow instruments.

When to use it: every wrapped load, handpieces and hollow instruments — anywhere a gravity sterilizer cannot guarantee steam penetration.

KEY FEATURES & BENEFITS

- Pre-vacuum removes air so steam reaches every surface
- Post-vacuum drying — loads come out dry, ready to store
- 134°C and prion cycles on a clear LCD display
- Four-button operation — staff trained in minutes
- Health Canada & ETL/CSA approved

Monitor with ATLAS Bowie-Dick, Type 5 and rapid BIs

Ask us about financing terms, installation and bundle pricing with ATLAS monitoring products.



A proud Canadian company · Health Canada registered products
info@penguinhealthgroup.com · +1 905-901-1579 · penguinhealth.ca

CE CREDITS AVAILABLE

Specialized expertise. Practical support.

A compliant sterilizer does not guarantee a compliant practice. Our IPAC visits find the gaps — before an inspector, an insurer or a patient does.

FREE
IPAC COURSES FOR
ATLAS SUMMIT
CLIENTS

COURSE 1

3.5-Hour In-Depth Visit

A full review of your sterilization room and team, with a written report you can file for compliance.

- 90 min on-site assessment
- 60 min team education
- 30 min hands-on coaching at the point of work
- 30 min review of findings
- Written compliance report

COURSE 2 · MOST POPULAR

2-Hour Focused Training

Targets one specific concern — a new hire, a new protocol or a recent finding.

- 1 hour assessment of the concern
- 1 hour education for the team members involved
- Delivered in your own sterilization room

What we review

Six areas, in your own sterilization room — not a classroom.

Workflow & layout

How instruments move from dirty to sterile to storage.

Sterilizer monitoring

Biological, chemical and Bowie-Dick practice.

Cleaning & packaging

Ultrasonic, washer and pouch technique.

Failures & procedures

What happens when a test fails, and whether the team knows it.

Documentation

What is recorded, where, and whether it is retrievable.

Team readiness

Whether every operator can explain the "why," not just the "how."

Course pricing

3.5-HOUR ASSESSMENT & TRAINING

~~\$3,500~~ **\$1,500** with Summit Bundle

2-HOUR READINESS ASSESSMENT

~~\$2,000~~ **FREE** with Summit



ATLAS™

Elevating your sterilization to new heights.

ATLAS RAPID

ATLAS SHIELD

ATLAS WASH

ATLAS TRACKER

A proud Canadian company · Health Canada registered products

info@penguinhealthgroup.com · +1 905-901-1579 · penguinhealth.ca

Penguin Health Incorporated · Mississauga, Ontario, Canada · **Lifetime warranty on all our products**

Specifications are subject to change. Product availability and permitted claims vary by market. ATLAS is a trademark of Penguin Health Incorporated.



LUNA
WHITENING

NEW • LUNA PRACTICE STARTER BUNDLE

3 Take Home Kits + 1 Desensitizing Gel + 1 Gingival Barrier

Everything to start offering whitening this week

REG \$491.95

SAVE \$92.95

\$399

LUNA WHITENING TAKE HOME KIT



A NEW ADVANCED TAKE HOME TOOTH WHITENING SYSTEM WITH LED LIGHT TECHNOLOGY

INCLUDED

- 2 WHITENING SYRINGES X 5ML
- ACTIVATION LIGHT
- SILICONE TRAY
- GUIDE & INSTRUCTIONS



3+1

NET \$98.99

3+1

REG \$131.99 /KIT

REORDER CODE
PH-LUNATHK

LUNA WHITENING IN OFFICE FULL KIT



SINGLE PATIENT KIT - 30% CP + 30% HP

OUR STRONGEST GEL YET!

- ONE 40 MINUTE SESSION PATIENTS SEE A CHANGE OF 2-10+ SHADES WITH NO CUSTOM TRAYS.
- INCLUDED 3ML OF PRE-WHITENING PROPHY PASTE + 3ML OF POST CARE WHITENING GEL.



4+2

NET \$87.99 /KIT

4+2

REG \$131.99 /PACK

REORDER CODE
PH-MAXX

LUNA DESENSITIZING GEL (5ML SYRINGE)

- STRENGTHENS AND REMINERALIZES TOOTH ENAMEL, REDUCING SENSITIVITY AFTER WHITENING.
- PROMOTES ACTIVE REMINERALIZATION WITH FLUORIDES, CALCIUM, AND PHOSPHATE.
- EASY TO APPLY, EFFECTIVE FOR NIGHTTIME USE.
- CONTAINS SODIUM FLUORIDE TO BIND MINERALS TO ENAMEL.



3+1

NET \$26.99 /SYRINGE

3+1

REG \$35.99 /SYRINGE

REORDER CODE
PH-DG5S

LUNA WHITENING IN OFFICE PACK



2 WHITENING SYRINGES X 3ML - 30% CP/HP



3+1

NET \$35.99 /KIT

3+1

REG \$47.99 /PACK

REORDER CODE
PH-LTH1

LUNA

GINGIVAL (GUM) BARRIER SYRINGE
(3ML - TREATS 3 CLIENTS)



2+1

NET \$39.99 /SYRINGE

2+1

REG \$59.99 /SYRINGE

REORDER CODE
PH-LGBS3

PROFESSIONAL ADVICE

Dental Professionals - Select the right LUNA Whitening Concentration by assessing the patient's oral health and sensitivity. Use a lower concentration (16% CP to 25% CP) for sensitive teeth and higher concentrations (35% HP) for faster results. Tailored treatments ensure safe, effective, and satisfiable whitening! (Devila et al, 2020)*



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



PRECISION EDGE SYSTEMS SPECIALS

PENGUIN PRECISION EDGE KIT SECTIONAL MATRIX SYSTEM KIT

CONTAINS:

100 MATRICES - 25 EACH
3.5 mm - 4.5 mm - 5.5 mm - 6.5 mm

75 WEDGES - 25 EACH

Small - Medium - Large

- 1 UNIVERSAL RING
- 1 NARROW RING
- 1 FORCEPS & 1 PIN TWEEZERS

50%
OFF



NET \$349.50

INTRO KIT A

50% OFF

REG. \$699 /KIT

REORDER CODE

PH-EDGEKIT

PENGUIN PRECISION EDGE RING REFILL

COMPARES TO PALODENT RING
MADE WITH NITI MATERIALS

- 1 UNIVERSAL RING | PH-DMB0495
- 1 NARROW RING | PH-DMB0496

MORE FLEXIBLE AND
DURABLE THAN
COMPETITOR!

1 PER BOX

2+1



NET \$83.32 /1BX

2 + 1

REG. \$124.99 /BX

PENGUIN PRECISION TULWAR MATRICE WEDGE REFILL

- CAN BE DOUBLE-WEDGED
- THE TIP ANGLED TO PROTECT RUBBER DAM
- EASY TO INSERT AND REMOVE

SMALL/DARK BLUE | PH-DIW0903
MEDIUM/LIGHT BLUE | PH-DIW0902
LARGE/WHITE | PH-DIW0901

COMPARES TO PALODENT
MATRICE WEDGE

BOX OF 100

2+1



Small



Medium



Large

NET \$39.99 / BX OF 100

2 + 1

REG. \$59.99 /BX

PENGUIN PRECISION MATRICE REFILL

- SUBGINGIVAL MATRICES AVAILABLE!

COMPARES TO PALODENT
MATRICE

3.5MM | PH-DMB0811
PRE-MOLAR - 4.5MM | PH-DMB0812
MOLAR - 5.5MM | PH-DMB0813
MOLAR W/EXTENSIONS - 6.5MM | PH-DMB0814

BOX OF 50



3.5mm



4.5mm



5.5mm



6.5mm

NET \$46.66 / BX OF 50

2 + 1

REG. \$69.99 /BX



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



O-STAR 7 Modes

Unprecedented
Curing Light
incorporating
various strengths

Charging base with
built-in light meter

- 7 MODES (NORMAL, HIGH, TURBO, ORTHO, SOFT, PULSE, AND CARIES CHECKER)
- INDUSTRY-LEADING PERMEABILITY AND LIGHT-LOSS REDUCTION!
- SLIM-HEAD AND 360-DEGREE ROTATION
- WIDE-RANGE WAVELENGTH CURING!

O-STAR CORDLESS CURING LIGHT

Q4 SPECIAL PRICE

\$959

each • PH-WP-O-STAR

REGULAR PRICE

\$1,999

YOU SAVE

\$1,040

REORDER CODE

PH-WP-O-STAR

BUY WITH A HARVARD BUNDLE

EXTRA 25% OFF

O-Star for just **\$719.25**

WWW.PENGUINHEALTH.CA

KEY SELLING POINTS • AT A GLANCE

Harvard — German quality since 1892



MultiChrome Universal

PH01855 compules • PH01854 syringe

- Chameleon effect for a perfect shade match
- Spherical micro-filler technology
- Compules or 3 g syringe

NET
\$151.99
4+1
REG
\$189.99
20 compules ×
0.25 g


MultiChrome Flow

PH-7100003

- Multi-shade flowable — blends into the tooth
- Pairs with MultiChrome Universal
- Made in Germany

NET
\$111.99
4+1
REG
\$139.99
2 × 1 mL
syringes


UltraFill

Shade codes on p.24

- Permanent high gloss
- Opalescent, tooth-like fluorescence
- Low polymerization shrinkage

NET
\$119.96
4+1
REG
\$149.95
Syringe pack

BioCal-CAP

Codes on p.26

- Bioactive resin-modified MTA cement
- pH 11 promotes hydroxyapatite formation
- Light-cure: continue treatment immediately

NET
\$87.99
4+1
REG
\$109.99
1 g syringe

InterLock ONE

Codes on p.27

- Universal: works with every etching technique
- Only one thin layer needed
- Optional self-cure activator

NET
\$187.49
3+1
REG
\$249.00
5 mL bottle

ZirconCore

PH-7083599 (A2) • PH-7083600 (A3)

- Nano-zirconia reinforced
- Core build-up AND post cementation
- Grinds like dentin, without grooves

NET
\$142.49
3+1
REG
\$189.99
5 mL pack
+25% OFF

Buy any Harvard bundle and get an **extra 25% off the O-Star curing light** — just \$719.25 (page 19).



HARVARD®

Harvard MultiChrome

- Advanced micro filler technology, spherical fillers
- Color pigments
- Chameleon effect for perfect match with substance
- Flexible logistics

NEW



Cuts like dentin

Made in
Germany



EN



CLICK OR SCAN QR CODE FOR
BROCHURES & RESEARCH!



★ KEY FEATURES

- ✓ Chameleon effect: perfect shade match
- ✓ Spherical micro-filler technology
- ✓ Made in Germany

MULTICHROME UNIVERSAL COMPULES

20 COMPULES X
0.25G/PACK!



REORDER CODE

PH01855

NET \$151.99 /PK

4+1

REG \$189.99 /PACK

MULTICHROME UNIVERSAL 3G SYRINGE



REORDER CODE

PH01854

NET \$119.20 /PK

4+1

REG \$149.00 /PACK

MULTICHROME BLOCK N MASK COMPULES

BUY MULTICHROME
AND GET SPECIAL PRICE!

BOX OF 20 COMPULES



REORDER CODE

PH-7100007

NET \$151.99 /PK

4+1

REG \$189.99 /PACK

MULTICHROME UNIVERSAL FLOW SYRINGE PACK 2X1ML



REORDER CODE

PH-7100003

NET \$111.99 /PK

4+1

REG \$139.99 /PACK



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca






A PROUD CANADIAN
COMPANY

Harvard UltraFill

Ultra Gloss, Opalescent, High Performance

Composite for anterior and posterior restorations.



UF Harvard UltraFill
 HM Heliomolar
 DU DuraFill VS
 GR Gradia Anterior
 GAe G-aenial Anterior
 ES Esthet-X HD
 TE Tetric EvoCeram
 PF Harvard PremiumFill

Harvard

UltraFill

Ultra Gloss, Opalescent, High Performance Composite for **anterior** and **posterior** restorations

- Permanent high gloss
- Opalescence
- Tooth-like fluorescence
- Easy to work with
- Superior handling properties
- Low polymerization shrinkage
- High abrasion resistance
- Optimal physical properties
- Radiopaque

HEALTH CANADA
LICENSE
MDALL: 109285

1892
seal
HARVAR



IntroKit S, 3 x 3 g



3 g



20 x 0.25 g

HARVARD ULTRAFILL COMPULES

NET \$151.96 /PK

4+1

REG \$189.95 /PACK


HARVARD ULTRAFILL SYRINGE SPECIALS

NET \$119.96 /PK

4+1

REG \$149.95 /PACK

Universal Opacity easy technique



Universal



CLICK OR SCAN QR CODE FOR
BROCHURES & RESEARCH!

- Selected shades in **one** opacity
- In 80 % of clinical cases best compromise of opacity

Harvard UltraFill „Universal Opacity“	3 g syringe
Shade	
A1 U	7111111
A2 U	7111112
A3 U	7111113
A3.5 U	7111114
BL U	7111151
Harvard Applier OptiTips®	

★ KEY FEATURES

- ✓ Permanent high gloss
- ✓ Opalescent, tooth-like fluorescence
- ✓ Low polymerization shrinkage
- ✓ Universal opacity: best match in 80% of cases



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



Harvard UltraFill Flow

Ultra Gloss, Opalescent, Flowable Composite
for **anterior** and posterior restorations

- >> Easy polishability, permanent high gloss
- >> Opalescence
- >> Optimal flow consistency
- >> Thixotropic
- >> Tooth-like fluorescence
- >> Reduced polymerization shrinkage
- >> High abrasion resistance
- >> Optimal physical properties
- >> Easy to work with
- >> Radiopaque
- >> Available in ComforTip®F

HEALTH CANADA
LICENSE
MDALL: 109285



Order Details

2 x 1 ml syringe		20 x 0.25 g ComforTip® F	
A1	7112111	A1	7112211
A2	7112112	A2	7112212
A3	7112113	A3	7112213
A3.5	7112114		
B1	7112121		
BL	7112151		
Incl. 8 needle tips H18			
Harvard NeedleTips H18		7095158	
Refill bag with 50 needle tips			
Harvard Applier OptiTips®		7095200	

No stock item, please contact the Harvard sales team.



20 x 0.25 g

NET \$119.99_K


4+1

REG. \$149.99 /PACK

NET \$151.99_{/PK}

4+1

REG. \$189.99 /PACK



ComforTip®F 0.25 g

★ KEY FEATURES

- ✓ Thixotropic — stays where placed
- ✓ Easy to polish, permanent high gloss
- ✓ Bendable 40° ComforTip

Dentist MSc. Andreas Kluschke



1: Tooth 14 and 15 Wedge-shaped defects



2: Class V restorations after high gloss polish

„Easy and quick to polish, durable high gloss.“



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



Harvard BioCal®-CAP

Bioactive, light cure, resin-modified MTA cement
for direct and indirect pulp capping

- >> For direct and indirect pulp capping in the treatment of deciduous and adult teeth
- >> Calcium release and a high pH value (pH 11) promote the formation of hydroxyapatite as well as tertiary dentine
- >> Creates an environment hostile to bacteria
- >> Moisture tolerant
- >> Virtually no solubility
- >> Good Radiopacity
- >> Fast, after light cure treatment can be continued immediately
- >> Thixotropic properties

HEALTH CANADA
LICENSE
MDALL: 109285

ORDER DETAIL

1 g Syringe, 12 Needle Tips

7081550

NET \$87.99 /PK

4+1

REG. **\$109.99** /PACK

PACK OF 1 SYR

NET \$178.40 /PK


4+1

REG. **\$223** /PACK


PACK OF 4 SYR



Direct pulp capping



Indirect pulp capping



Significant release of bioactive calcium promotes the formation of hydroxyapatite and tertiary dentine.



The plaster for the pulp.



★ KEY FEATURES

- ✓ pH 11 drives hydroxyapatite formation
- ✓ Light-cure: continue immediately
- ✓ Direct & indirect pulp capping



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



Harvard InterLock® ONE

Universal Adhesive

HEALTH CANADA
LICENSE
MDALL: 109285



- >> Universal: suitable for all tooth conditioning techniques: with or without etching gel, wet or dry
- >> Only one thin layer necessary
- >> High and reproducible bond strength to enamel and dentine even under long term stress
- >> Effective marginal seal for low risk of postoperative sensitivity
- >> Excellent performance also with self or dual cure composites
- >> Convenient and precise application with the triangle ergonomic bottle with defined drop size

If light curing of Harvard InterLock® ONE cannot be achieved, please mix with Harvard InterLock® ONE Self Cure Activator.

Harvard InterLock® ONE Self Cure Activator

- >> Activator for safe adhesion in indications, where light curing of Harvard InterLock® ONE cannot be ensured (e.g. in the root canal)
- >> Ideal for cementing of posts with dual and self cure composites

NET \$187.49

3+1

REG \$249



5 ml



50 x 0.1 ml

NET \$300.79

4+1

REG \$375.99

Order Details	
Universal Adhesive	
5 ml bottle	7083611
50 x 0.1 ml SingleDose	7083613
Self Cure Activator	
5 ml bottle	7083612
Microbrush®	
Refill bag with 50 applicators	7095156



5 ml

NET \$172.49

3+1

REG \$229.99



Harvard InterLock®

Very strong, self-etch, light cure bonding within two steps

- >> Long-term, strong and safe adhesion to enamel and dentine
- >> High and even bond strength
- >> Extremely reliable
- >> Easy and fast application
- >> Excellent seal
- >> Not technique and moisture sensitive
- >> Hydrophilic for wet bonding technique
- >> For perfect margins to protect against secondary caries

The gold standard

★ KEY FEATURES

- ✓ Every etching technique
- ✓ One thin layer only

NET \$319.20

4+1

REG \$399



2 x 5 ml

Order Details	
2 x 5 ml bottle	
50 Microbrush®, 5 Appli-Pads	7083603
Microbrush®	
Refill bag with 50 applicators	7095156



HEALTH CANADA
LICENSE
MDALL: 109285

Harvard ZirconCore

Dual cure core build-up and post cementation composite



- Excellent consistencies for both indications
core build-up and **post cementation**
- Very good mechanical properties for permanent restorations
- Contains a substantial amount of Nano-zirconia
- Fluoride release
- Thixotropic properties
- Very good radiopacity
- Self-mixing material with the Minimix syringe
- For optimal results, use with
Harvard InterLock® ONE /
Harvard InterLock® ONE Self Cure Activator
or with **Harvard Bond SE Dual**

Grinding without grooves.

★ KEY FEATURES

- ✓ Nano-zirconia reinforced
- ✓ Core build-up AND post cementation
- ✓ Grinds like dentin, no grooves

Excellent
mechanical
properties
in all
curing
modes.

3+1

NET \$142.49

3+1

REG \$189.99

REORDER CODE - A2

PH-7083599

REORDER CODE - A3

PH-7083600

ORDER INFORMATION

SHADE A2 - 5ML PACK**SHADE A3 - 5ML PACK**

orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca

A PROUD CANADIAN
COMPANY

GS - STERILIZATION POUCHES

CLASS 4 SELF-SEALING

200/BOX - DUAL INDICATOR


PROMO 7+3

BRAND	#SIZE	REG.	NET.
GS POUCHES	2.75" X 10"	\$12.99	\$9.09
GS POUCHES	3.5" X 10"	\$12.99	\$9.09
GS POUCHES	3.5" X 6.25"	\$11.99	\$8.39
GS POUCHES	5.25" X 11"	\$24.99	\$17.49

CHEETA

BARRIER FILM ROLL

CLEAR & BLUE



NET \$24.74/ROLL

PROMO 3+1

REG \$32.99/ROLL

CLEAR

REORDER CODE

PH-BF02C

OPTIM OS1 WIPE

1 MINUTE KILL TIME

160/CAN



REORDER CODE

PH-OSICA-W12

INTRODUCING NEW HALYARD WRAPS

SMS STERILIZATION WRAP



SIZE	BOX	PRICE	#REORDERS	10% DISCOUNT
12X12	500	\$79.95	PH-D31212	\$71.95
15X15	500	\$116.49	PH-D31515	\$104.84
20X20	500	\$168.95	PH-D32020	\$152.05
24X24	500	\$224.99	PH-D32424	\$202.49

INTRODUCING HB CSR WRAP



SIZE	BOX	PRICE	#REORDERS	10% DISCOUNT
15X15	500	\$89.99	PH-HBCSR1515	\$80.99
20X20	500	\$149.99	PH-HBCSR2020	\$134.99
24X24	500	\$198.99	PH-HBCSR2424	\$179.09



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



A PROUD CANADIAN
COMPANY

OPTIM® 33 TB SURFACE WIPES

CAN OF 160 WIPES

BUY 12

BUY 24



REORDER CODE

PH-SCI33W12

BUY 24

NET\$16.49 /CAN OF 160 WIPES

BUY 12

NET\$18.95 /CAN OF 160 WIPES

OPTIM 33TB LIQUID 4L

BUY 4



NET \$49.49 /BTL

BUY 4

REG \$59.99L/BTL

REORDER CODE

PH-SCI334X4

MARK3 GENERAL PURPOSE

ULTRASONIC SOLUTION 4L

BUY 4



NET \$39.99 /JUG

BUY 4

REG \$49.99 /JUG

REORDER CODE

PH-100-0092

SURFACE DISINFECTANT GERMINOL RTU

4L BOTTLE



BUY 4



NET \$39.99 /BTL

BUY 4

REG \$59.99 /1 BTL

REORDER CODE

PH-GERRTU-4L

BIOENZYMATIC ULTRASONIC CLEANER

83036 64 TABLETS/BOX

5+1



NET \$58.32 /BOX

5+1

REG \$69.99 /BOX

REORDER CODE

M3100-7640

PENGUIN LIQUID HANDSOAP

4L BOTTLE



30%
DISCOUNT



NET \$27.99 /BTL

30% DISCOUNT

REG \$39.99 /BTL

REORDER CODE

PH-PHS-4L



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



A PROUD CANADIAN
COMPANY

INTCO ● NITRILE GLOVES

S M L XL

BOX OF 100



NET \$6.49/BOX
10+10
REG \$12.99/BOX OF 100

S	REORDER CODE SY-100S	REORDER CODE SY-100M	M
L	REORDER CODE SY-100L	REORDER CODE SY-100XL	XL

INTCO BLACK ● NITRILE GLOVES

S M L XL

BOX OF 100



NET \$8.99/BOX
7+3
REG \$12.85/BOX OF 100

S	REORDER CODE PH-ADV100SB	REORDER CODE PH-ADV100MB	M
L	REORDER CODE PH-ADV100LB	REORDER CODE PH-ADV100XLB	XL

PRI-MED ● BLUEBERRY NITRILE GLOVES

S M L XL

BOX OF 250



NET \$20.99/BOX
7+3
REG \$29.99/BOX OF 200

S	REORDER CODE PM6-2806	REORDER CODE PM6-2807	M
L	REORDER CODE PM6-2808	REORDER CODE PM6-2809	XL

RONCO SILKTEX LATEX EXAM GLOVES - 5MIL

XS S M L

BOX OF 100



NET \$13.99/BOX
7+3
REG \$19.99/BOX OF 100

XS	REORDER CODE PH-MT-819XS	REORDER CODE SK-819S	S
M	REORDER CODE SK-819M	REORDER CODE SK-819L	L

ZIRCONIA CUTTING DIAMOND BUR CROWN & ENDO ACCESS

\$89.99/PACK OF 5

NET \$71.99 / PK



ZB6LLC (Z856L/018C)



ZB2MLC (Z882/014C)



ZALL3SC (Z801L/014C)



ZC12SSNC (Z830/010C)



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



A PROUD CANADIAN
COMPANY

BRONX

AIR WATER SYRINGE TIPS

PACK OF 250

4+2

NET \$32.66/PACK

4+2

REG \$48.99 /PK OF 250

REORDER CODE

BRO-2381



ZEN

SALIVA EJECTORS
CLEAR WITH WHITE TIP
CLEAR WITH BLUE TIP

PACK OF 100

8+2

NET \$4.76/PACK

8+2

REG \$5.95 /PK OF 100

WHITE TIP

REORDER CODE

ZEN-2384

BLUE TIP

REORDER CODE

ZENCB100



CURVED UTILITY

#412 SYRINGE 12CC

PACK OF 50

BUY 4+

NET \$29.99/BAG

BUY 4+

REG \$39.99 /PK OF 50

REORDER CODE

100-5212



HB DENTAL BIBS

BLUE

PACK OF 500

3+1

NET \$37.49/CASE


3+1

REG \$49.99 /PK OF 250

PURPLE

REORDER CODE

PH00003



HB TRAY COVER

WHITE

CASE OF 1000

3+1

NET \$37.49 /CASE

3+1

REG \$49.99 /CASE OF 1000

REORDER CODE

HBWTP-1000



HB TIPS

MIXING TIPS

BLUE - 1:1 - REORDER CODE - 100-1418

PINK - 5.0MM - 100-1417 - REORDER CODE - 100-1417

YELLOW - 4.2MM - REORDER CODE - M3100-1412

PURPLE - 7.5MM - REORDER CODE - M3100-1416

M3 HIGH PERFORMANCE MIXING TIPS

(GREEN - 0.5MM - HEAVY BODY & BITE REGISTRATION)

- REORDER CODE - PH-100-1413

NET \$29.69/BOX

10% DISCOUNT

REG \$32.99/BOX

10% DISCOUNT



CHEETA PROPHY CUPS

LATCH TYPE

CHE-PCLS-BL SOFT - BLUE
CHE-PCLS-PI SOFT - PINK

PACK OF 144

3+1

NET \$27.74/PACK

3+1

REG \$36.99/PACK

REORDER CODE

PH00019

REORDER CODE

PH-CHE-PCLS-PI



SCREW TYPE

CHE-PCS-BL SOFT - BLUE
CHE-PCS-PI SOFT - PINK

PACK OF 144

3+1

NET \$24.74/PACK

3+1

REG \$32.99/PACK

REORDER CODE

PH-CHE-PCS-BL

REORDER CODE

PH00020



SNAP-ON

CHE-PCSOS-B SOFT - BLUE

PACK OF 144

3+1

NET \$14.99/PACK

3+1

REG \$19.99/PACK

REORDER CODE

PH-CHE-PCSOS-B



MILO

PRE-BENT APPLICATOR TIPS

18 GAUGE - PINK | REORDER | MIL-2387

20 GAUGE - BLACK | REORDER | MIL-2388

22 GAUGE - GREY | REORDER | MIL-2389

25 GAUGE - BLUE | REORDER | MIL-2390

4+1

NET \$11.19/BAG

4+1

REG \$13.99/BAG OF 100



CHEETA HEADREST COVER

CASE OF 500 - 10X13

3+1

NET \$40.49/CASE

3+1

REG \$53.99/CASE

BLUE

REORDER CODE

PH00011

WHITE

REORDER CODE

PH00010



KAI MICRO APPLICATORS

REGULAR | REORDER | KAI-2304

FINE | REORDER | KAI-2302

SUPERFINE | REORDER | KAI-2301

BOX OF 400

3+1

NET \$26.21/BOX

3+1

REG \$34.95/BOX OF 400
4 TUBES OF 100



MARK 3

INTRAORAL MIXING TIPS

INTRAORAL TIPS YELLOW | REORDER | 100-0060

INTRAORAL TIPS CLEAR | REORDER | 100-0059

20% DISCOUNT

PACK OF 100

NET \$27.99/PACK

20% DISCOUNT

REG \$34.99/PACK OF 100



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



A PROUD CANADIAN COMPANY

PT-A™

AIR POLISHER & ULTRASONIC SCALER

Elevate your practice today!



CLICK OR SCAN THIS QR CODE FOR VIDEOS AND RESOURCES!

0%

Financing For 24 months

\$369.98/MONTH

*terms and conditions apply

PROMO

GET 20 POWDER BOTTLES FREE

REORDER CODE

PH01992



PROPHYLAXIS POWDERS

FAST

Sodium Bicarbonate
(Supra - 65um)



GENTLE MINI

Glycine
(Sub & Supra - 25um)



Recommended angle:
30°-60°



Recommended distance:
3-5mm



MOBILE CART

SPECIAL

\$989/UNIT

REORDER CODE

PH-WP-MC



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



A PROUD CANADIAN COMPANY

WOODPECKER PT-A / PT-B

WOODPECKER PT-A / PT-B
SUPRAGINGIVAL HANDPIECE
(AP-1) - AP-1

3+1

NET \$449.96 /UNIT
(MIX AND MATCH)

3+1

REG \$599.95 /UNIT



REORDER CODE

PH-WP-AP-1

WOODPECKER PT-A / PT-B
SUBGINGIVAL HANDPIECE
(AP-2) - AP-2

3+1

NET \$449.96 /UNIT
(MIX AND MATCH)

3+1

REG \$599.95 /UNIT



REORDER CODE

PH-WP-AP-2

WOODPECKER PT-A / PT-B
ULTRASONIC HANDPIECE
(TZ-2L) - TZ-2L

NET \$449.96 /UNIT
(MIX AND MATCH)

3+1

REG \$599.95 /UNIT



REORDER CODE

PH-WP-TZ-2L

WOODPECKER PT-A / PT-B
SUBGINGIVAL AIR POLISHING
NOZZLES

2+1

NET \$46.63 /PACK

2+1

REG 69.95 /PK OF 10 NOZZLES



REORDER CODE

PH-WPNOZZLE10

WOODPECKER DTE | ULTRASONIC SCALING TIPS

NET \$69.95/TIP

P33 TIP



REORDER CODE

PH-P33

P50L TIP



REORDER CODE

PH-P50R

P50R TIP



REORDER CODE

PH-P59

P56 TIP



REORDER CODE

PH-P50L

P59 TIP



REORDER CODE

PH-P56

E60 TIP



REORDER CODE

PH-E60

2+2

NET \$69.95 /TIP

REG \$139.95 /TIP

WOODPECKER GENTLE

(GLYCINE - SUPRAGINGIVAL 65UM -
200G BOTTLE - PROPHYLAXIS POWDER)

NET \$33.32 /BTL

PH-WP-GENTLE

2+1

REG \$49.99 /120G BTL



WOODPECKER GENTLE MINI

(GLYCINE - SUBGINGIVAL + SUPRAGINGIVAL
POWDER 25UM - 120G BOTTLE -
PROPHYLAXIS POWDER)

NET \$33.32 /BTL

PH-WP-GENTLE-MINI

2+1

REG \$49.99 /120G BTL



WOODPECKER FAST

(GLYCINE - SUPRAGINGIVAL 65UM -
200G BOTTLE - PROPHYLAXIS POWDER)

NET \$33.32 /BTL

PH-WP-FAST

2+1

REG \$49.99 /120G BTL



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



PENGUIN

CLASSIC LOW SPEED HANDPIECE

1:1 Direct Drive	2+1	REG \$299.95 /UNIT
	NET \$199.97 /UNIT	2+1
(A1008) G100-S1H1 • For CA burs (ø2.35) • Latch type • Max Speed: 30,000 min		
REORDER CODE PH-A1008		

1:1 Direct Drive	2+1	REG \$299.95 /UNIT
	NET \$199.97 /UNIT	2+1
(A1009) G100-S1H2 • For CA burs (ø2.35) • Clean Head System • Push Button Chuck • Max Speed: 40,000 min		
REORDER CODE G100-S1H2		

	2+1	REG \$299.95 /UNIT
NET \$199.97 /UNIT	2+1	
(A1012) G100-ST2 • For HP burs (ø2.35) • Stainless Steel Body • Clean Head System • Max Speed: 40,000 min		
REORDER CODE PH-G100-ST2		

PROPHY HANDPIECE

4:1 Reduction	1+1	REG \$399.95 /UNIT
	NET \$199.99 /BOX	1+1
(A1017) G200-SC • For screw-in cups and brushes • Stainless Steel Body • 360° Rotation • w/o Spray • Paste Defense System • Max Speed: 5,000 min		
REORDER CODE PH-G200-SC		

4:1 Reduction	1+1	REG \$399.95 /UNIT
	NET \$199.99 /BOX	1+1
G200-LA • For CA burs (ø2.35) • Push Button Chuck • Stainless Steel Body • Clean Head System • Ball Bearing • Max Speed: 5,000 min		
REORDER CODE PH-G200-LA		

SLOW SPEED AIR MOTOR



2+1

REG **\$499.95**/UNIT

NET \$333.30/UNIT



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



ORAL SURGERY INSTRUMENTS



Free instrument must be
of same or lesser value

SCALERS ONLY
BUY 9+3



HYGIENE & EXAM INSTRUMENTS

BUY 12 +3

Free instrument must be
of same or lesser value

BUY 5

GET 1

FREE

PRODUCT NAME	ITEM CODE	PRICE
204S Sickle Scaler	1184	34.99
H6/H7 Hygienist Scaler	1175	34.99
Super Taha Scaler	STAHA	34.99
Taha 1 Anterior Sickle Scaler	TAHA1	34.99
Taha 2 Posterior Sickle Scaler	TAHA2	34.99
Taha 3 Posterior Sickle Scaler	TAHA3	34.99
Taha 4 Posterior Sickle Scaler	TAHA4	34.99
Super Taha Scaler	STAHA	34.99
11/12 Gracey Curette	1234	34.99
13/14 Gracey Curette	1235	34.99
15/16 Gracey Curette	1236	34.99
1/2 Gracey Curette	1229	34.99
1/2 Mini Rigid Gracey Curette	G12MR	34.99
1/2 Arhum Gracey Curette	AG12	34.99
11/14 Arhum Gracey Curette	AG1114	34.99
13/14 Rigid Gracey Curette	G1314R	34.99
15/16 Rigid Gracey Curette	G1516R	34.99
7/8 Younger Good Universal Curette	1212	34.99
4R/4L Columbia Universal Curette	1216	34.99
13/14 Columbia Universal Curette	1217	34.99
H5/J33 Hygienist Jacquette Sickle Scaler	1176	34.99
5/6 Langer Universal Curette	1206	34.99
1/2 Langer Universal Curette	1204	34.99
13/14 Indiana University Universal Curette	1208	34.99
ODU 11/12 Elongated Explorer	1122AF	29.99
ODU 11/12 Explorer	1122	29.99
Signature Mirror Handle, Cone Socket	1168	21.99
College Utility Dressing Plier	1157	19.99
Meriam Utility Dressing Plier	1158	19.99
23/UNC15 Expro	1151	29.99
23/QW Williams Expro	1146	29.99
23/CP 12 Siyah Lachak Expro	1150S	39.99
UNC15/11.5B Probe	1137	29.99
18 Ball Burnisher, 1.6MM 2.4MM	1918	32.99
1 Goldstein Plastic Filling Instrument, Siyah Series	2467	42.99
2 Goldstein Plastic Filling Instrument, Siyah Series	2468	42.99
3 Goldstein Extra Flex Plastic Filling Instrument, Siyah Series	2469	42.99
4 Goldstein Extra Flex Plastic Filling Instrument, Siyah Series	2470	42.99
27/29 Ball Burnisher, 1.6MM 3.1MM, Siyah Series	300-104	42.99
Small/Medium Ball Burnisher, Siyah Series	300-132	44.99
T-Flex Composite Instrument	TFLEX	98.99
E2 Spoon, 1.5MM Excavator	1808	32.99
#18 Spoon, 1.5MM Excavator	1813	32.99
26/27S Ball Burnisher, 1.3MM 2.1MM	1919	32.99
28/29 Ball Burnisher	1921	32.99
L7 Levy Carver	2397	32.99
1/2 Hollenback Carver, Siyah Series	1904S	42.99
Interproximal Carver, Long, Siyah Series	2369	42.99

HiTeck

Your Partner In Dentistry

FLIGHT A12 OPERATORY CHAIR PACKAGE

(INCLUDES FREE DENTAL STOOL)

\$599/month

0%

FINANCING FOR
48 MONTHS



Features:

- A12 Hydraulic Patient Chair
- State of the Art Hands free Sensor Controlled Flight Torch LED Light
- Synchronized Chair and Back Movement for superior Ergonomics
- Choose from 9 Plush Sewn or Seamless Upholstery Colors and 2 Ultraleather Colors
- Self Contained Water Bottle System
- Thin Backrest and Streamlined Seat rest
- Designed and Engineered in Canada
- DCI Compatible Parts

AUTO CLAVE B

\$330/month

The **Clave B** Autoclave features a pre/post vacuum that provides full penetration of steam into wrapped, porous, and hollow loads. The pre-vacuum phase enables the complete removal of air from the chamber while the post-vacuum air removal ensures that the load is completely dry.

The combination of the pre/post vacuum allows closed door drying and sterilization for the most difficult loads.

The Clave B is an intuitive and easy to use Class B Sterilizer. All operations are controlled with only 4 buttons, with 3 buttons that control the program selection and 1 to start the sterilization cycle.

Health Canada and ETL/CSA Approved
Excellent Sterilizer for use in Dental Offices, Veterinarian Offices, Medical Offices, Podiatry, Chiropodist, Tattoo and Beauty Salons

0%

FINANCING FOR
36 MONTHS



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



★ Q4 SPECIAL PRICE ★

REG \$19,999

\$14,999

YOU SAVE \$5,000

Smile 2 Intraoral Scanner · PH-SMILE2

PLUS — INCLUDED FREE

✓ High-spec laptop **\$2,000 value**

✓ Scanning software **\$3,000 value**

✓ 100 GB cloud storage **FREE**

TOTAL PACKAGE VALUE \$24,999+

YOU SAVE OVER \$10,000



SMILE2 Intraoral Scanner


SMART CLOUD, EFFORTLESS ACCESS

Utilize a high-capacity cloud system for the upload and secure storage of data. This system ensures seamless integration with laboratories, facilitating the one-click transmission of scan data. Advanced partner management tools further optimize collaboration and enhance operational efficiency. Cases can be shared instantly through QR codes, providing patients with convenient access to digital models. This functionality significantly improves communication between dental professionals and their patients.




Data Hub: Simplified Management

Utilize a comprehensive data statistics system to track business performance with clear charts, ensuring you maintain control over your business's growth.



Versatile Platforms Availability Plus Offline Capability

Support for PC, web browsers, and mobile devices enhances flexibility. Support offline functionality to scan and access data anytime, anywhere.



PROUD TO BE CANADIAN

BUY LOCAL SUPPORT LOCAL



**PENGUIN
HEALTH**

A Proud Canadian Company

YOUR TURNKEY OFFICE PARTNERS

CONTACT US TODAY

OUR PROCESS

- Building Solutions
- Project Management
- Value - Added Approach

CONSTRUCTION PROMO

FREE Consultation, **FREE** Office
Design, and **Discount**
on Equipment!



orders@penguinhealthgroup.com



905-901-1579



www.penguinhealth.ca



Notes

orders@penguinhealthgroup.com • 905-901-1579 • penguinhealth.ca



ATLAS™

Elevating Your Sterilization
to New Heights.

ONE SYSTEM • ONE CANADIAN SUPPLIER

Every link in sterilization, *elevated.*



ATLAS Rapid
20-minute results



ATLAS Shield BI
Ultra rapid steam



ATLAS Shield Type 5
Integrating indicator



ATLAS Wash-U-Test
Cleaning verification



ATLAS Tracker
Digital records

IN-CLINIC EXCLUSIVE

Book your ATLAS Summit Bundle visit

Save over \$5,000, including a FREE IPAC in-clinic course.
Your representative brings the bundle to you.

ORDER & QUESTIONS

905-901-1579

orders@penguinhealthgroup.com
penguinhealth.ca